

#### **Open Access** Review



# Harnessing monocyte-derived dendritic cells for evaluating T cell response by mixed leukocyte reaction

Honorio Torres-Aguilar<sup>\*</sup><sup>®</sup>, Sorely Adelina Sosa-Luis<sup>®</sup>, Alexia Almaraz-Arreortua<sup>®</sup>, William de Jesús Ríos-Ríos<sup>®</sup>

Clinical and Basic Immunology Research Department Biochemical Sciences Faculty, Universidad Autónoma "Benito Juárez" de Oaxaca, Oaxaca City 68120, Mexico

\*Correspondence: Honorio Torres-Aguilar, Clinical and Basic Immunology Research Department Biochemical Sciences Faculty, Universidad Autónoma "Benito Juárez" de Oaxaca, Oaxaca City 68120, Mexico. qbhonorio@hotmail.com Academic Editor: Jinming Han, Capital Medical University, China Received: December 4, 2024 Accepted: June 5, 2025 Published: June 30, 2025

**Cite this article:** Torres-Aguilar H, Sosa-Luis SA, Almaraz-Arreortua A, Ríos-Ríos WdJ. Harnessing monocyte-derived dendritic cells for evaluating T cell response by mixed leukocyte reaction. Explor Immunol. 2025;5:1003201. https://doi.org/10.37349/ei.2025.1003201

## Abstract

The mixed leukocyte reaction (MLR) is a pivotal in vitro assay for evaluating T-cell responses stimulated by allogeneic antigen-presenting cells (APCs). Dendritic cells (DCs) are the most efficient stimulatory cells. However, the scarcity of circulating DCs in peripheral blood limits their isolation for research or clinical use. In contrast, monocytes, which are abundant and easily accessible, can be differentiated into monocytederived DCs (moDCs) in vitro and have emerged as the most practical and efficient stimulatory cells for MLR due to their accessibility and robust allostimulatory capabilities. This review aims to describe the scientific rationale and evidence for using moDCs in MLR assays to assess T-cell alloreactivity. Its methodology outlines the protocols for experimental, preclinical, and biosafety assays that have demonstrated the practicality of moDCs in evaluating and quantifying the alloresponse of naïve and memory CD4+ and CD8+ T cells, as well as the effects of immunomodulatory factors, immune monitoring, and tolerogenic strategies in the context of transplantation. Additionally, it illustrates how moDC-mediated MLRs have provided critical insights into understanding alloimmunity processes and antigen-specific T-cell responses in cancer immunotherapy, autoimmune diseases, and vaccine development, with potential implications for personalized medicine and immunotherapy optimization. In conclusion, despite ongoing challenges such as standardization and scalability in massive cell production, the current understanding and reproducible results of moDC applications in MLRs highlight their potential to develop innovative strategies focused on immune monitoring.

## **Keywords**

Monocyte-derived dendritic cells, mixed leukocyte reaction, alloreactivity, antigen-specific responses, cancer, autoimmunity, vaccines, immunotherapy

© The Author(s) 2025. This is an Open Access article licensed under a Creative Commons Attribution 4.0 International License (https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, sharing, adaptation, distribution and reproduction in any medium or format, for any purpose, even commercially, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made.



# Introduction

Under steady-state conditions, dendritic cells (DCs) take up commensal microbes and apoptotic cells to maintain immunologic tolerance. Conversely, when DCs are exposed to an inflammatory environment, this exposure triggers DC activation, transforming them into specialized cells that process, present, and activate pathogen-derived antigen-specific T cells [1]. This process strengthens the immune response and facilitates the induction of antibody-producing plasma cells [2].

Surface molecule expression and single-cell RNA sequencing (scRNA-seq) have characterized various DC subsets in human blood and other tissues [3]. Nevertheless, blood-circulating DCs are sparse, exhibiting a semi-mature state or migrating toward their final location and differentiation [4]. Therefore, isolating them from peripheral blood may yield insufficient numbers, rendering them ineffective as stimulatory antigen-presenting cells (APCs) for in vitro assays. Fortunately, the bloodstream contains monocytes that can differentiate into DCs [3]. In vivo, monocytes serve as a continuous reserve of precursors for tissue-resident macrophages and monocyte-derived DCs (moDCs) [5]. Since the method for obtaining in vitro-produced moDC was described in 1994 by Sallusto and Lanzavecchia [6], numerous variations of this method have emerged, demonstrating successful moDC production for both research and clinical applications.

As efficient APCs, in vitro-produced moDCs have been evaluated in preclinical strategies where the immune system requires antigen-specific stimulation [e.g., tumor-specific antigens (TSA) in cancer] [7] or inactivation (e.g., alloantigens before transplantation) [8]. Given their accessibility and efficiency as stimulatory cells, moDCs have been utilized in mixed leukocyte reaction (MLR) assays, which represent an effective in vitro strategy for quantifying T cell allo- and antigen-specific responses, as well as the effects of immunomodulatory agents [9–15].

## moDCs as inducers of T cell alloresponse

The success and efficiency of an MLR in evaluating the potential T cell response against specific alloantigens depend on selecting the appropriate APCs. Blood is the most readily available source of human mononuclear cells, including monocytes and B cells [16]. However, while these cells can effectively present alloantigens, their numbers vary among individuals and may not represent the most suitable APCs involved in graft rejection [17, 18]. DCs are specialized APCs known to be the most efficient cells at inducing T cell alloresponse; however, their numbers in peripheral blood are limited [4]. This section analyzes the characteristics of moDCs that render them the most suitable and efficient APC for MLR.

### DCs: modulating the fate of the immune system

The immune system's functional complexity depends on maintaining homeostasis with the commensal microbiota, regulating opportunistic microbes, preventing pathogenic microorganisms, and aiding in tissue regeneration in the event of injury [19]. To this end, in a steady state, tissue-resident immature DCs exhibit high phagocytic capacity, facilitating the uptake of commensal microbes and apoptotic self-cells resulting from normal tissue repair. Under these conditions, DCs show low expression of costimulatory molecules (CD80, CD86, and CD40), inflammatory and migratory chemokine receptors such as CC chemokine receptor 7 (CCR7), major histocompatibility complex (MHC) class II molecules, and pro-inflammatory cytokines, while also displaying constitutive expression of these DCs toward secondary lymphoid tissues to present innocuous microbe-derived or self-antigens to T cells induces an unresponsive state (anergy) or generates regulatory T cells (Tregs), maintaining immune tolerance and avoiding potentially harmful chronic inflammation [23].

On the other hand, an inflammatory milieu caused by tissue damage provokes the maturation of DCs due to pathogen-associated molecular patterns (PAMPs) present in microorganisms and damage-associated molecular patterns (DAMPs) derived from injured cells [20]. This alarming environment activates metabolic, cellular, and genetic transcription signaling pathways that, beyond a microbicidal function,

enable DCs to specialize in processing pathogen-derived proteins into antigenic peptides and loading them onto MHC molecules [24]. In addition, due to an overexpression of CCR7 induced by inflammation, DCs accelerate their migration toward secondary lymphoid tissues [21], where changes in cytokines and costimulatory profiles result in antigen presentation and the recognition by specific T cells that produce effector lymphocytes, thereby strengthening the immune response [25] and working together to stimulate antibody-producing plasma cells [26].

#### **Human DC subsets**

The expression of cytokines, along with the profiles of costimulatory and inhibitory molecules, varies among resident and migrating DCs across different tissues [27]. Similarly, each DC subset responds differently to pro-inflammatory and anti-inflammatory stimuli, reflecting its unique functionality and heterogeneity [28]. Due to the limited availability of human tissues, researchers have mainly characterized the DC subpopulations using mouse models [29]. In humans, peripheral blood DCs and certain tissue-resident DCs from experimentally accessible organs, such as the skin, tonsils, and spleen, are well characterized. However, our understanding of most human DC subsets has evolved from phenotypic, molecular, and functional correlations with mouse DCs [30].

Various human DC classifications have been established based on the expression of surface molecules. Additionally, scRNA-seq technology, which assesses the transcriptomic and functional profiles, has enabled the description of six DC groups in human peripheral blood: DC1 to DC6. Each subset is endowed with specific transcription factors and markers [3]. However, the traditional classification defines two main groups in human peripheral blood: conventional DCs (cDCs), which include cDC1 and cDC2, and plasmacytoid DCs (pDCs). All of these can migrate to different tissues. Interestingly, in peripheral blood, there are also monocytes capable of differentiating into DCs [21].

#### In vivo human moDCs

In peripheral blood, monocytes are a heterogeneous group of immune cells derived from the myeloid lineage in the bone marrow, transiently circulating in the bloodstream, and representing about 2–10% of the total white blood cell count. These partially differentiated cells have short half-lives in blood circulation, approximately one day. Currently, three subsets of blood monocytes have been identified in humans based on the differential expression of CD14 (Toll-like receptor 4 coreceptor) and CD16 (Fc $\gamma$  receptor III). CD14++/CD16- classical monocytes exhibit high phagocytic abilities, making up 80–90% of total monocytes. In contrast, CD14++/CD16+ or intermediate monocytes and CD14+/CD16++ or non-classical monocytes together account for around 10% of total monocytes. Due to their location in the marginating pool within blood vessels and their primary role in tissue remodeling, they are referred to as "patrolling" monocytes [31].

Given their high numbers and functional properties, classical monocytes serve as a rapid and continuous reserve of precursors for tissue-resident macrophage populations during homeostasis and inflammation, and potentially, for moDCs [32]. However, exploring moDC in human tissues has posed a significant challenge due to ethical and logistical difficulties. It has become clear that, in vivo, most DCs originate from a specific blood-circulating precursor called pre-DC, which is distinct from monocytes [33]. For instance, the study of patients with mutations in the *IRF8* and *GATA2* genes revealed a deficiency in blood monocyte numbers and functions, a lack of certain dermal DC subtypes, and a reduced number of macrophages, while the numbers of Langerhans cells (epidermal DCs) remained intact [34]. This evidence indicates that, under homeostatic conditions, certain dermal DCs are derived directly from circulating monocytes (steady-state moDCs). On the other hand, under inflammatory conditions, cells recruited to the skin due to monocyte extravasation maintain phenotypic characteristics of DCs and lack the macrophage markers CD16 and CD163. They have received various designations such as moDCs, inflammatory DCs, monocyte-like inflammatory cells (Infl mo-like), inflammatory macrophages (Inf mac), or monocyte-derived macrophages (Monomac) [35]. Additionally, strong evidence from animal models, human transplantation, marker analysis, and disease contexts shows that monocytes differentiate into DCs and

macrophages in vivo, particularly during inflammation. While some tissue-resident macrophages originate from embryonic sources, monocytes serve as crucial responders in dynamic or pathological conditions [3, 15, 36].

#### In vitro generation of moDCs

In addition to describing the significant role of DCs in modulating the immune system, experimental strategies have emerged to harness their potential for clinical applications [7]. However, being aware of the limited numbers (0.1% to 1.0% of peripheral blood mononuclear cells) [37], the challenges of isolating blood-circulating DCs [38], and the impossibility of separating them from other human tissues [39], in vitro differentiated DCs derived from blood monocytes have served as the primary source of human DCs for research and cellular therapy [8, 40].

Sallusto and Lanzavecchia first described the methodology for differentiating monocytes into cells with the morphological, phenotypic, and functional features of DCs in 1994 by culturing monocytes in the presence of granulocyte-macrophage colony-stimulating factor (GM-CSF) and interleukin (IL)-4 [6]. Since then, numerous alternatives to DC precursors, purification methods, culture mediums, differentiation factors, cytokine concentrations, maturation stimuli, and phenotypic characterization have been successfully developed [41, 42].

The primary challenge in generating moDCs for biomedical research or applying them to cellular therapy is establishing and validating viable production procedures that meet good manufacturing practice requirements to produce functional APCs capable of priming quantifiable T-cell-mediated responses [9]. The consensus reveals that obtaining antigen-loaded mature moDC may take a time-consuming protocol of 7–9 days (at least five days to differentiate immature moDC with high phagocytic abilities, one day to allow uptake of the antigen and processing, and followed by two additional days to obtain mature moDC). Although some proposals suggest reducing the entire production time to 48–36 hours with encouraging experimental evidence, they still leave the expectation of obtaining reproducible results and clinical applications [43].

#### **Current clinical applications of moDCs**

The immunotherapy concept encompasses clinical strategies focused on modulating the immune system. Under this premise, DCs are the most efficient immune cells when antigen-specific activation (e.g., against TSAs in cancer or pathogen-specific antigens in infection) or inactivation (e.g., against autoantigens in autoimmune disease (AD) or alloantigens before transplantation) is desired to induce either immune overactivation or immunosuppression, respectively [44].

As described above, the in vitro production of moDCs allows for the harnessing of their functional plasticity by using pro-inflammatory or anti-inflammatory factors during their differentiation and maturation to obtain immunogenic (Imm-moDCs) or tolerogenic (Tol-moDCs) moDCs [7]. In this regard, Imm-moDCs have been used in interventional clinical trials for conditions such as melanoma, bladder cancer, primary advanced carcinoma of the oral cavity, squamous cell carcinoma of the head and neck, and HIV infection, among others (https://clinicaltrials.gov/). On the other hand, Tol-moDCs have been evaluated to assess the effects of intervention in multiple sclerosis, rheumatoid arthritis, musculoskeletal diseases, type 1 diabetes, and joint disease (https://clinicaltrials.gov/), as well as kidney transplant rejection [40]. Nevertheless, despite the vast amount of information, the first and only U.S. Food and Drug Administration-approved moDC-based cellular therapy is a vaccine, sipuleucel-T, or Provenge, which has been in use since 2010 [45].

#### moDC involvement during in vivo alloresponse

Monocytes are continuously produced in the bone marrow, where they differentiate into DCs or macrophages in various tissues. As a result, they constantly migrate toward injured tissues, drawn by inflammatory chemokines (CCL2, CCL7, CCL5, CCL20, CX3CL1) across the endothelium. Their migration is influenced by maturation stimuli, including TNF- $\alpha$ , IL-1, microbial products, and CD40L on activated T cells,

platelets, and mast cells. Additionally, they can take up both self- and pathogen-derived antigens [5, 36]. However, while monocytes have high phagocytic abilities to capture dead cells, they primarily degrade ingested proteins instead of correctly processing them to produce appropriate peptides for antigen presentation. Nevertheless, their differentiation into moDCs enhances processing capacities for presenting proteins derived from dead cells [46].

Therefore, although performed under sterile conditions, transplantation itself—due to the inflammatory environment produced by tissue damage—causes the maturation and migration of both the donor (embedded in the allograft) and the recipient (in adjacent tissues) DCs, including moDCs attracted by the chemokines and induced by the differentiation factors produced in the inflammatory milieu [47]. Hence, T cells from the graft recipient may recognize allogeneic MHC molecules via two pathways mediated by DCs. A) Through a so-called "direct pathway" involving the recognition of unprocessed MHC molecules in the graft. Donor DCs are instigated to migrate toward the recipient's secondary lymphoid tissues. Because positive T cell selection restricts them to recognizing only self-MHC molecules, there is a high frequency of alloantigen-specific recipient naïve and memory T cells that recognize allogeneic MHC molecules expressed on donor DCs [48]. B) The "indirect pathway," although described as less significant for graft rejection, involves the entry of recipient DCs into the grafted organ. After taking up and processing donor MHC and other alloantigens from dead cells, these DCs migrate toward recipient secondary lymphoid tissues to activate alloantigen-specific recipient cytotoxic T cells, which attack the graft, along with helper T cells that collaborate to induce alloantibody-producing plasma cells [49, 50]. moDCs are also implicated in chronic rejection, resulting in fibrosis and the development of chronic allograft injury in liver transplantation through their interactions with T cells and other immune cells [51].

#### moDCs application to reduce in vivo alloresponse

Conventional immunosuppressive therapies focus on affecting the functionality of T and B cells [52]. However, counteracting the immune responses that lead to graft rejection becomes a significant challenge once these cells have been activated [53]. Therefore, it is anticipated that DC depletion in the graft and DC inactivation in the host may impede alloantigen presentation [18], and therapeutic strategies to modify DC functions to harness their tolerogenic properties would enhance the efficiency of organ transplantation [40, 54]. However, considering their inter-tissue location, depleting DCs in the grafted organ is not feasible, and depleting DCs in the host could lead to immunodeficiency [55]. Although specific strategies to modulate DC activation have been developed, their ineffectiveness in exclusively affecting in vivo DCs limits their clinical application [56, 57]. Therefore, using DCs to prevent graft rejection is regarded as cell therapy through the production of autologous Tol-moDCs [58].

Tol-moDCs can be produced in vitro using protocols that require immunomodulatory agents and genetic technology, such as anti-inflammatory cytokines (IL-10, TGF- $\beta$ ) [59], approved immunosuppressive drugs [60], experimental inhibitory molecules (vitamin D3) [61], genetic engineering (antisense oligonucleotides targeting costimulatory molecules or promoting anti-inflammatory factors) [62–64], or a combination of several strategies. Yet, no consensus exists on the best procedure to produce therapeutic-grade Tol moDCs [65, 66].

The hallmark of Tol moDCs is low expression of costimulatory molecules (mainly CD40) and proinflammatory factors (IL-12p70) [67]. Tol-moDC can induce T cell tolerance by promoting anergy, enhancing T regs, facilitating apoptosis, or skewing phenotypes through cell-cell contact-dependent mechanisms or soluble factors [8, 23, 40].

Currently, autologous Tol-moDCs have been used in phase I/II monocentric clinical trials to evaluate biosafety for patients with renal insufficiency receiving a kidney transplant (ClinicalTrials.gov ID: NCT06243289 and NCT02252055), aiming to induce transplant tolerance. Further clinical trials are needed to assess doses, frequency, and administration routes, among other clinical parameters. The preclinical assays in rhesus monkeys indicate that the tolerogenic effect of Tol moDC is mediated by the induction of Treg [68].

# MLR application to evaluate alloreactive memory T-cell responses in the context of transplantation

Before a transplant, a patient must undergo several tests to maximize the safety of the procedure. These evaluations include blood group and HLA typing, as well as the search for unexpected allospecific antibodies. Currently, increasing graft compatibility primarily relies on HLA matching performed through molecular methods [69]. Nevertheless, despite those theory-matching tests, the existence of donor alloantigen-specific T cells in the receptor signifies an enormous challenge in transplantation. The above is because given the T cell restriction to recognize just self-MHC molecules—even in the absence of previous events of allogenic stimulation (transfusion, transplants, or pregnancies)—each person possesses a high frequency of alloreactive naïve and, even worse, memory helper (CD4+) and cytotoxic (CD8+) T cells with the potential to become activated and instigate graft rejection [17, 70, 71]. Since the MLR serves as an in vitro assay to analyze effective T-cell responses primed for alloantigen recognition [72–74]. Therefore, in addition to the theory matching tests mentioned, MLR may be conducted to assess the functional histocompatibility of the donor-recipient pair.

Several biosafety and preclinical trials have demonstrated the feasibility of moDCs for evaluating and quantifying the alloresponse of preexisting memory T cells prior to a transplant, along with the impact of inhibitory factors, such as suppressive drugs, on T cell responses (ClinicalTrials.gov ID: NCT03726307, NCT06243289, NCT00495755, NCT00213707).

Applying the MLR to evaluate the alloresponse involves assessing the likelihood of graft rejection by coculturing recipient-derived T cells, purified from peripheral blood, as responders stimulated with donor-derived moDCs [11]. On the other hand, the procedure may involve coculturing peripheral blood-purified donor-derived T cells primed with recipient-derived moDCs to analyze the possibility of graft-versus-host disease (GvHD) [75–78] (Figure 1).

MLR can be performed using peripheral blood mononuclear cells (PBMC) purified by differential centrifugation with Ficoll from the whole blood of the donor and recipient [79]. However, under these conditions, APC function is performed by monocytes and B cells, whose numbers may vary among patients. These cells do not represent the most suitable APCs involved in graft rejection and for MLR with effective results. A high density of costimulatory and MHC class II molecules is crucial to triggering an efficient MLR. Although monocytes and B cells possess these molecules, their expression levels are low, which provides them with weak allostimulatory capabilities [16].

Depending on the specific transplanted tissue, responder cells for an MLR may be obtained from different sources, albeit with certain limitations. For example, conducting a one-way MLR with bone marrow-derived responder cells indicates that, within this tissue, most cells are not T cells, and many others demonstrate persistent basal proliferation [80]. On the other hand, purifying lymphocytes from the thymus, lymph nodes, or spleen as responder cells is appropriate for MLR in animal models [81, 82]. However, it is not feasible in human assays, and stimulator cells must be obtained from PBMC or induced in vitro [10-12, 83].

The MLR can be two-way when both T cell populations, from the donor and recipient, proliferate [84] or one-way, when only one of those lymphocyte populations can respond [85]. Although the first provides a general view of the donor-recipient's histocompatibility, it does not allow for the evaluation of each response. Therefore, the latter is preferable because it enables the evaluation of the likelihood of graft rejection by quantifying the proliferation of only recipient-derived T cells as responders or the possibility of GvHD by analyzing donor-derived T cells. This is achieved by abolishing the proliferative ability of stimulator cells through irradiation [86] or treating them with mitomycin C [85].

Some recipients of stem cells or kidney allografts—because of their inherent condition—are often immunocompromised, and their peripheral blood T cells may be affected. Therefore, in this case, to evaluate whether these cells have the potential to proliferate, the inclusion of HLA-DR mismatched controls should be analyzed. Then, the response in the MLR between donor and recipient is compared with the reaction toward the health controls [83].



Figure 1. Methods using moDC for evaluating T cell responses in MLR. (1) PBMC are obtained by centrifugation gradients, such as Ficoll, from peripheral blood or leukapheresis products. (2) In the MACS system, monocyte and T cell populations are purified by positive selection using anti-CD14-, anti-CD4-, or anti-CD8-conjugated magnetic microbeads. (3) Monocytes are differentiated in the presence of GM-CSF and IL-4 for 4 to 7 days to obtain imoDCs. (4) mmoDCs are obtained by inducing their maturation with PAMPs (e.g., LPS, Poly I: C) or inflammatory stimuli (e.g., LPS, TNF, IL-6, PGE<sub>2</sub>) for an additional 2 days. (5) The efficiency of the obtained moDC can be analyzed by flow cytometry, which includes differentiation (CD1a, CD83) and maturation (CD80, CD86, CD40) markers. (6) An MLR can evaluate the likelihood of graft rejection through coculturing recipientderived T cells stimulated with donor-derived moDCs. (7) A second MLR may analyze the possibility of graft-versus-host disease, incorporating donor-derived T cells primed with recipient-derived moDCs. (8) In the MLRs, T cells recognize alloantigens processed and presented on allogeneic MHC molecules by allogeneic APCs. The highest quantifiable T cell proliferation is observed around days 5-8, with proliferation and activation state analyzed by flow cytometry using the CFSEdilution method and measuring CD25 and CD69 expression, respectively. Additionally, analyzing soluble factors (cytokines and chemokines) in supernatants provides information about the alloresponse intensity. moDC: monocyte-derived dendritic cells; MLR: mixed leukocyte reaction; PBMC: peripheral blood mononuclear cells; GM-CSF: granulocyte-macrophage colonystimulating factor; imoDCs: immature moDCs; mmoDCs: mature monocyte-derived dendritic cells; PAMPs: pathogen-associated molecular patterns; PGE<sub>2</sub>: prostaglandin E<sub>2</sub>; MHC: major histocompatibility complex; APCs: antigen-presenting cells; CFSE: carboxy fluorescein succinimidyl ester

MLR application in immunosuppressive drug testing involves using moDCs in preclinical studies to evaluate the in vitro effect of immunosuppressive drugs on alloreactivity, such as calcineurin and mTOR inhibitors, as well as co-stimulatory blockade agents on both DC stimulatory abilities and T cell activation, proliferation, and cytokine profiles [51, 87–91]. On the other hand, in addition to methods such as flow cytometry for analyzing T cell activation, monitoring donor-specific antibodies, and using specific radiotracers to visualize T cell trafficking, among others [92], moDC can be used in immune monitoring in clinical settings to track specific T cell responses to alloantigens during the post-transplant phase. This procedure assesses the frequency and intensity of T-cell alloresponse, providing insights into the risk of rejection or tolerance [93].

## moDCs application to evaluate T cell response in cancer immunotherapy

moDCs are widely used to assess the immune responses triggered by cancer immunotherapies, such as allogeneic tumor cell-based vaccines or adoptive tumor-specific T-cell transfer. In these clinical interventions, moDCs are utilized in vivo to present TSA to T cells and subsequently in vitro through an MLR to monitor the resulting immune activation toward TSA, allo, and autoantigens [9, 10, 94–100].

The MLR can evaluate the strength of tumor-specific T-cell responses. This is accomplished by mixing T cells from the patient with moDCs previously loaded with TSA or tumor-derived lysates. The resulting cellular events—proliferation, activation, and cytokine profiling—provide information about tumor immunogenicity and the immune system's ability to recognize and trigger a tumor-specific immune response or even the risk of autoimmunity [100].

Adoptive cell therapy (ACT) involves the infusion of autologous, tumor-specific T-cells that have been expanded ex vivo or genetically engineered to target cancer cells. In this immunotherapy, MLR is used to measure whether these transferred cells recognize, proliferate, and kill tumor cells [101]. However, in allogeneic ACT, MLR may also predict the risk of unwanted immune alloresponse, such as donor T cell rejection by the recipient's immune system or GvHD mediated by donor T cells that damage recipient cells. This could compromise the therapy's success, which depends on the compatibility and persistence of donor T cells that must recognize and destroy host tumor cells [102].

Additionally, including moDCs loaded with different TSAs in various MLR assays allows for the identification of the most immunogenic TSA that elicits strong T-cell responses as potential candidates for developing targeted immunotherapies, such as personalized cancer vaccines or ACT [98–100]. Likewise, the likely alterations of T-cell activation mediated by the tumor microenvironment (TME) can be evaluated by adding tumor-derived factors (e.g., soluble components or immune cells from the tumor) to the MLR, providing insights into how immunotherapies could counteract these effects [98, 100]. Finally, the MLR can be developed for personalized immunotherapy approaches tailored to the specific immune reactivity of the tumor, utilizing autologous moDCs pulsed with their tumor-derived antigens [94].

# moDCs-based MLR applications to assess the immune profiling in autoimmunity

ADs are characterized by the pathological activation of preexisting low-affinity T cells that recognize selfantigens. Depending on their maturation state, DCs have been involved as both instigators and protectors of AD [103]. Therefore, as the best APC, moDCs-based MLR provides a suitable tool to investigate the presence and contributions of autoreactive T cells for the diagnosis, treatment, and monitoring of these diseases (e.g., rheumatoid arthritis, multiple sclerosis, type 1 diabetes, systemic lupus erythematosus, etc.).

Self-antigen-loaded moDCs (e.g., myelin in multiple sclerosis, GAD65 in type 1 diabetes,  $\beta$ 2GPI in antiphospholipid syndrome, or joint proteins in rheumatoid arthritis) are utilized in MLR to quantify the presence of autoreactive T cells and analyze their activation state and cytokine profile (e.g., IFN- $\gamma$ , TNF- $\alpha$ , IL-6, IL-17) in autoimmune patients, as well as assess the functionality of Tregs in maintaining immune tolerance [59, 103–105].

moDCs-based MLR enables the distinction of each specific T cell polarization contribution (e.g., Th1, Th17, or Treg cells) to autoimmune pathology. MLR can also be employed to examine the Th1/Th17 dominance. In many ADs, IL-17 (Th17 cells) or IFN- $\gamma$  (Th1 cells) are critical cytokines driving chronic inflammation and tissue damage [106].

MLR can also be applied to analyze the role of individual pro- and anti-inflammatory cytokines, environmental factors (e.g., toxins), pathogen-derived antigens, immunosuppressive drugs (e.g., corticosteroids, methotrexate), biologic agents (e.g., TNF inhibitors), drugs targeting specific immune checkpoints (such as CTLA-4 or PD-1), tolerance-inducing vaccines, or peptide-based therapies in the beginning, progression, and counteraction of the autoimmune process [107] and simultaneously investigate

the signaling pathways activated in T cells, such as the NF-kB, JAK/STAT, and MAPK pathways, which are frequently dysregulated in autoimmune conditions [108]. Therefore, since the immune response can differ among AD patients, the MLR can be used to predict how a specific AD patient might respond to a particular treatment in order to create tailored therapeutic strategies for that individual [109, 110].

# moDC-based MLR in vaccine development

The MLR applications in vaccine development encompass evaluating antigen immunogenicity, biosafety, and personalized vaccine design. This allows for the simultaneous assessment of T cells and cytokine profiles to optimize vaccine formulations [94].

moDC-based MLR can be used to analyze T cell activation and proliferation, as well as polarization [CD4+ helper T cells (Th1, Th2, Th17), CD8+ cytotoxic T cells, Tregs], and cytokine secretion (e.g., IFN- $\gamma$ , IL-2, TNF- $\alpha$ ) against specific candidate antigens for a vaccine. Given this, the immune balance will determine whether this antigen may induce an immune response with cell memory and antibody production, while avoiding excessive immune tolerance or even autoimmunity [111].

Likewise, in combination with vaccine antigens, MLR may be applied to optimize vaccine formulations, test different adjuvants, evaluate various vaccine delivery systems such as liposomes, nanoparticles, or viral vectors, assess cross-reactivity (i.e., whether T cells might respond to similar antigens), identify the most immunodominant epitopes within a vaccine antigen, and measure the recall response of memory T cells to subsequent doses [112–114].

Since the MLR involves human moDC-T cell interactions, it represents a better humanized preclinical model than animal models. This assay may provide more accurate predictions of how a vaccine would perform in humans. The MLR can be adapted to assess immune responses under different genetic backgrounds or in individuals with immunocompromised conditions (e.g., elderly subjects or patients with cancer), chronic diseases (e.g., diabetes, obesity), co-infections (HIV), or even specific T cell responses in children. These variations ensure that vaccine formulations can be adjusted for different immune system profiles, enhancing safety and efficacy across all populations [114–116].

# Methods using moDC for evaluating T cell responses in MLR

A better approach to assessing a preexisting T-cell response should consider that several DC subsets, with specific functions and tissue distribution near the inflammatory process, may contribute differentially to antigen presentation [18]. However, due to their accessibility and efficiency, only in vitro-produced moDCs have been extensively used as the best APCs (stimulatory cells) in MLR assays [9–12, 83].

The most common source for obtaining monocytes is peripheral blood or leukapheresis products centrifugation gradients (Ficoll) separate PBMC. Monocytes are purified by magnetic cell sorting (MACS), which achieves high purity rates above 95%. Monocytes are differentiated in vitro by culturing them with GM-CSF and IL-4 for 4 to 7 days, generating immature moDCs (imoDCs). Then, maturation is induced by adding stimuli such as LPS, TNF, Poly I: C, IL-6, and PGE<sub>2</sub>, either alone or in different combinations, for 2 more days to obtain mature moDCs (mmoDCs). The status of the obtained cells is assessed by measuring the expression balance of differentiation/maturation and inhibitory markers such as CD1a, CD83, CD86, CD40, Tim-3, PD-L1, and HLA-DR using flow cytometry.

Responder cells can be represented by whole PBMC or partially purified naïve or CD3+ T cells. However, the alloresponse can be analyzed in greater detail by purifying CD4+, CD8+, and CD45RA-memory T cells using MACS.

The moDC: T cell cocultures are performed in a 1:10 ratio for 4 to 7 days. The balance between cell proliferation and death can be measured using simple methods such as cell counting and viability assessments through light microscopy and trypan blue staining in the Neubauer chamber, although these may yield indefinite and unquantifiable results [117]. Cell proliferation can be quantified by integrating radioactive tritium [<sup>3</sup>H]-labeled thymidine [118] into the MLR culture. Still, this method does not evaluate

cell death or activation and requires a scintillation counter. Therefore, flow cytometry represents an advanced technology for cell analysis and data acquisition, enabling the simultaneous analysis and quantification of cell proliferation [e.g., CFSE-dilution method, 5-bromo-2'-deoxyuridine (BrdU), Ki-67 expression, WST-8/CCK8 tetrazolium assay] related to CD4- or CD8-specific responses, cytokine expression (e.g., IFN-γ, IL-2), cell death (e.g., annexin V, propidium iodide, 7AAD), activation (e.g., CD25, CD69), and the induction of Tregs (e.g., Foxp3, CTLA-4) [13, 59, 104] and coculture supernatants can be utilized to quantify cytokines and chemokines related to the alloresponse intensity and T cell migratory abilities, such as IFN-γ, IL-12, CCL2, CCL3, CXCL9, and CXCL10 (Figure 1 and Table 1).

Responder cells	Source of moDCs	Days of stimulation	Proliferation method	Application	Reference and year
Peripheral blood CD4+	Peripheral blood	5 days	CD25 and Ki-67	To evaluate phenotypic changes associated with lymphocyte responsiveness	[9]
and CD8+ 1 cells purified by MACS	CD14+ monocytes purified by MACS		expression		2024
Peripheral blood naïve T cells purified by MACS	Peripheral blood CD14+ monocytes purified by MACS	7 days	CFSE dilution	To evaluate the effect of Sulfavant A in an MLR	[10]
					2023
PBMC	Monocytes from leukapheresis products	5 days	CFSE dilution	To evaluate recipient T cell responses facing donor- derived moDC before liver transplantation	[11]
					2021
Peripheral blood CD4+ and CD8+ T cells purified by MACS	Peripheral blood CD14+ monocytes purified by MACS	5 days	WST-8/CCK8 tetrazolium assay	To evaluate the effect of IL-35 on moDCs differentiation and stimulatory capability	[12]
					2018
Peripheral blood CD4+CD45RA– T cells purified by MACS	Peripheral blood CD14+ monocytes purified by MACS	5 days	CFSE dilution	To analyze tolerance induction on effector/memory T cells	[13]
					2010
Peripheral blood CD2+ T cells purified by MACS	Peripheral blood CD14+ monocytes purified by MACS	4 days	[ <sup>3</sup> H]-thymidine	To assess the	[119]
				impact of stably immature, donor-derived DC on alloimmune reactivity	2007
CD4+ peripheral blood T cells purified by MACS	Peripheral blood CD14+ monocytes purified by MACS	7 days	CFSE dilution and [ <sup>3</sup> H]- thymidine	To evaluate alloreactive T-cell responses before and after depletion of alloantigen- specific	[14]
					2004
				T cells	
CD3+ cells obtained from the peripheral blood	Peripheral blood CD14+ monocytes purified by MACS	5 days	[ <sup>3</sup> H]-methyl- thymidine	To monitor immunological and clinical response to DC vaccination	[15]
					2004

Table 1. Methods to anal	vze T cell responses	using moDCs as stimu	latory cells in an MI R
	yzo i oon rooponooo	aoning mobiol do otima	atory concentration

moDCs: monocyte-derived DCs; MLR: mixed leukocyte reaction; MACS: magnetic cell sorting; CFSE: carboxy fluorescein succinimidyl ester; <sup>3</sup>H: tritium; DC: dendritic cells; PBMC: peripheral blood mononuclear cells

In the MLR, T cells recognize alloantigens processed and presented on allogeneic MHC molecules by allogeneic APCs. This cell interaction (immunological synapse) induces that, among the entire T cell repertoire, only alloantigen-specific naïve and memory T cells become activated and proliferate. The highest measurable cellular proliferation occurs around 5–8 days [9, 10, 12, 83].

# Conclusions

Given their feasibility for in vitro generation and high allostimulatory capability, moDCs are the most efficient stimulatory cells in a mixed lymphocyte reaction. These properties support their application in evaluating preexisting memory T-cell alloresponse to analyze the likelihood of transplant success and the effects of immunomodulatory factors, immune monitoring, and tolerogenic strategies. Furthermore, the use of moDC in MLRs has provided critical insights into alloimmunity processes relevant to cancer immunotherapy, ADs, and vaccine development, facilitating the optimization of immunotherapeutic interventions.

# Abbreviations

AD: autoimmune disease APCs: antigen presenting cells CCR7: CC chemokine receptor 7 cDCs: conventional dendritic cells DCs: dendritic cells GM-CSF: granulocyte-macrophage colony-stimulating factor GvHD: graft-versus-host disease Imm-moDCs: immunogenic monocyte-derived dendritic cells MACS: magnetic cell sorting MHC: major histocompatibility complex MLR: mixed leukocyte reaction moDCs: monocyte-derived dendritic cells PAMPs: pathogen-associated molecular patterns PBMC: peripheral blood mononuclear cells scRNA-seq: single-cell RNA sequencing Tol-moDCs: tolerogenic monocyte-derived dendritic cells Tregs: regulatory T cells TSA: tumor-specific antigens

# **Declarations**

### Acknowledgments

The authors thank the Laboratorio Nacional de Citometría (LabNalCit), Oaxaca for their support for the moDCs-based MLR analysis from our laboratory, which has provided results for the writing of this review. During the preparation of this work, the authors utilized GitMind AI to create the hand elements of Figure 1. After using the tool, authors reviewed and edited the content as needed and will take full responsibility for the content of the publication.

### **Author contributions**

HTA: Conceptualization, Funding acquisition. SASL, AAA, WdJRR: Data curation, Investigation, Writing—original draft. All authors have read and approved the submitted version.

### **Conflicts of interest**

The authors declare that they have no conflicts of interest.

**Ethical approval** Not applicable.

**Consent to participate** Not applicable.

### **Consent to publication**

Not applicable.

#### Availability of data and materials

Not applicable.

### Funding

This research was funded by the Clinical and Basic Immunology Research Department Biochemical Sciences Faculty, Universidad Autónoma "Benito Juárez" de Oaxaca. A.A.A. has doctoral fellowships of CONAHCyT number 799779. The funders had no role in study design, data collection and analysis, decision to publish, or preparation of the manuscript.

## Copyright

© The Author(s) 2025.

## **Publisher's note**

Open Exploration maintains a neutral stance on jurisdictional claims in published institutional affiliations and maps. All opinions expressed in this article are the personal views of the author(s) and do not represent the stance of the editorial team or the publisher.

## References

- Bosma DMT, Busselaar J, Staal MD, Frijlink E, Mack M, Salerno F, et al. CD4<sup>+</sup> T-cell help delivery to monocyte-derived dendritic cells promotes effector differentiation of helper and cytotoxic T cells. Immunol Lett. 2025;275:107022. [DOI] [PubMed]
- Fayette J, Durand I, Bridon JM, Arpin C, Dubois B, Caux C, et al. Dendritic cells enhance the differentiation of naïve B cells into plasma cells in vitro. Scand J Immunol. 1998;48:563–70. [DOI] [PubMed]
- Villani A, Satija R, Reynolds G, Sarkizova S, Shekhar K, Fletcher J, et al. Single-cell RNA-seq reveals new types of human blood dendritic cells, monocytes, and progenitors. Science. 2017;356:eaah4573.
   [DOI] [PubMed] [PMC]
- 4. Orsini G, Legitimo A, Failli A, Massei F, Biver P, Consolini R. Enumeration of human peripheral blood dendritic cells throughout the life. Int Immunol. 2012;24:347–56. [DOI] [PubMed]
- Coillard A, Segura E. *In vivo* Differentiation of Human Monocytes. Front Immunol. 2019;10:1907.
  [DOI] [PubMed] [PMC]
- Sallusto F, Lanzavecchia A. Efficient presentation of soluble antigen by cultured human dendritic cells is maintained by granulocyte/macrophage colony-stimulating factor plus interleukin 4 and downregulated by tumor necrosis factor alpha. J Exp Med. 1994;179:1109–18. [DOI] [PubMed] [PMC]
- Plumas J. Harnessing dendritic cells for innovative therapeutic cancer vaccines. Curr Opin Oncol. 2022;34:161–8. [DOI] [PubMed]
- 8. Ochando J, Ordikhani F, Jordan S, Boros P, Thomson AW. Tolerogenic dendritic cells in organ transplantation. Transpl Int. 2020;33:113–27. [DOI] [PubMed] [PMC]
- 9. Mangelinck A, Dubuisson A, Becht E, Dromaint-Catesson S, Fasquel M, Provost N, et al. Characterization of CD4<sup>+</sup> and CD8<sup>+</sup> T cells responses in the mixed lymphocyte reaction by flow cytometry and single cell RNA sequencing. Front Immunol. 2024;14:1320481. [DOI] [PubMed] [PMC]
- Barra G, Gallo C, Carbone D, Ziaco M, Dell'Isola M, Affuso M, et al. The immunoregulatory effect of the TREM2-agonist Sulfavant A in human allogeneic mixed lymphocyte reaction. Front Immunol. 2023; 14:1050113. [DOI] [PubMed] [PMC]
- 11. Macedo C, Tran LM, Zahorchak AF, Dai H, Gu X, Ravichandran R, et al. Donor-derived regulatory dendritic cell infusion results in host cell cross-dressing and T cell subset changes in prospective living donor liver transplant recipients. Am J Transplant. 2021;21:2372–86. [DOI] [PubMed] [PMC]

- 12. Chen X, Hao S, Zhao Z, Liu J, Shao Q, Wang F, et al. Interleukin 35: Inhibitory regulator in monocytederived dendritic cell maturation and activation. Cytokine. 2018;108:43–52. [DOI] [PubMed]
- Torres-Aguilar H, Aguilar-Ruiz SR, González-Pérez G, Munguía R, Bajaña S, Meraz-Ríos MA, et al. Tolerogenic dendritic cells generated with different immunosuppressive cytokines induce antigenspecific anergy and regulatory properties in memory CD4+ T cells. J Immunol. 2010;184:1765–75.
   [DOI] [PubMed]
- Godfrey WR, Krampf MR, Taylor PA, Blazar BR. Ex vivo depletion of alloreactive cells based on CFSE dye dilution, activation antigen selection, and dendritic cell stimulation. Blood. 2004;103:1158–65.
  [DOI] [PubMed]
- Lee JJ, Kook H, Park MS, Nam JH, Choi BH, Song WH, et al. Immunotherapy using autologous monocyte-derived dendritic cells pulsed with leukemic cell lysates for acute myeloid leukemia relapse after autologous peripheral blood stem cell transplantation. J Clin Apher. 2004;19:66–70.
   [DOI] [PubMed]
- 16. Scheinecker C, Machold KP, Majdic O, Höcker P, Knapp W, Smolen JS. Initiation of the autologous mixed lymphocyte reaction requires the expression of costimulatory molecules B7-1 and B7-2 on human peripheral blood dendritic cells. J Immunol. 1998;161:3966–73. [PubMed]
- 17. Bingaman AW, Farber DL. Memory T cells in transplantation: generation, function, and potential role in rejection. Am J Transplant. 2004;4:846–52. [DOI] [PubMed]
- 18. Zhuang Q, Liu Q, Divito SJ, Zeng Q, Yatim KM, Hughes AD, et al. Graft-infiltrating host dendritic cells play a key role in organ transplant rejection. Nat Commun. 2016;7:12623. [DOI] [PubMed] [PMC]
- 19. Zoghi S, Masoumi F, Rezaei N. The immune system. In: Clinical Immunology. Elsevier; 2023. pp. 1–46. [DOI]
- 20. Quah BJC, O'Neill HC. Maturation of function in dendritic cells for tolerance and immunity. J Cell Mol Med. 2005;9:643–54. [DOI] [PubMed] [PMC]
- 21. Tiberio L, Prete AD, Schioppa T, Sozio F, Bosisio D, Sozzani S. Chemokine and chemotactic signals in dendritic cell migration. Cell Mol Immunol. 2018;15:346–52. [DOI] [PubMed] [PMC]
- 22. Carenza C, Calcaterra F, Oriolo F, Vito CD, Ubezio M, Porta MGD, et al. Costimulatory Molecules and Immune Checkpoints Are Differentially Expressed on Different Subsets of Dendritic Cells. Front Immunol. 2019;10:1325. [DOI] [PubMed] [PMC]
- Torres-Aguilar H, Blank M, Jara LJ, Shoenfeld Y. Tolerogenic dendritic cells in autoimmune diseases: crucial players in induction and prevention of autoimmunity. Autoimmun Rev. 2010;10:8–17. [DOI] [PubMed]
- 24. Leblanc-Hotte A, Audiger C, Chabot-Roy G, Lombard-Vadnais F, Delisle J, Peter Y, et al. Immature and mature bone marrow-derived dendritic cells exhibit distinct intracellular mechanical properties. Sci Rep. 2023;13:1967. [DOI] [PubMed] [PMC]
- 25. Bourque J, Hawiger D. Variegated Outcomes of T Cell Activation by Dendritic Cells in the Steady State. J Immunol. 2022;208:539–47. [DOI] [PubMed] [PMC]
- 26. Stebegg M, Bignon A, Hill DL, Silva-Cayetano A, Krueger C, Vanderleyden I, et al. Rejuvenating conventional dendritic cells and T follicular helper cell formation after vaccination. Elife. 2020;9: e52473. [DOI] [PubMed] [PMC]
- Segura E, Valladeau-Guilemond J, Donnadieu M, Sastre-Garau X, Soumelis V, Amigorena S. Characterization of resident and migratory dendritic cells in human lymph nodes. J Exp Med. 2012; 209:653–60. [DOI] [PubMed] [PMC]
- 28. Abboud G, Elshikha AS, Kanda N, Zeumer-Spataro L, Morel L. Contribution of Dendritic Cell Subsets to T Cell-Dependent Responses in Mice. J Immunol. 2022;208:1066–75. [DOI] [PubMed] [PMC]
- 29. Dutertre C. A new step in understanding mouse cDC ontogeny. Nat Immunol. 2024;25:383–4. [DOI] [PubMed]
- 30. Villar J, Segura E. Decoding the Heterogeneity of Human Dendritic Cell Subsets. Trends Immunol. 2020;41:1062–71. [DOI] [PubMed]

- 31. Ożańska A, Szymczak D, Rybka J. Pattern of human monocyte subpopulations in health and disease. Scand J Immunol. 2020;92:e12883. [DOI] [PubMed]
- 32. Ginhoux F, Jung S. Monocytes and macrophages: developmental pathways and tissue homeostasis. Nat Rev Immunol. 2014;14:392–404. [DOI] [PubMed]
- Capucha T, Mizraji G, Segev H, Blecher-Gonen R, Winter D, Khalaileh A, et al. Distinct Murine Mucosal Langerhans Cell Subsets Develop from Pre-dendritic Cells and Monocytes. Immunity. 2015;43: 369–81. [DOI] [PubMed]
- 34. Hambleton S, Salem S, Bustamante J, Bigley V, Boisson-Dupuis S, Azevedo J, et al. IRF8 mutations and human dendritic-cell immunodeficiency. N Engl J Med. 2011;365:127–38. [DOI] [PubMed] [PMC]
- 35. Mehta H, Angsana J, Bissonnette R, Muñoz-Elías EJ, Sarfati M. Inflammatory Skin Disorders: Monocyte-Derived Cells Take Center Stage. Front Immunol. 2021;12:691806. [DOI] [PubMed] [PMC]
- 36. Tang-Huau T, Segura E. Human in vivo-differentiated monocyte-derived dendritic cells. Semin Cell Dev Biol. 2019;86:44–9. [DOI] [PubMed]
- 37. Dzionek A, Fuchs A, Schmidt P, Cremer S, Zysk M, Miltenyi S, et al. BDCA-2, BDCA-3, and BDCA-4: three markers for distinct subsets of dendritic cells in human peripheral blood. J Immunol. 2000; 165:6037–46. [DOI] [PubMed]
- 38. Nair S, Archer GE, Tedder TF. Isolation and generation of human dendritic cells. Curr Protoc Immunol. 2012;99:7.32.1–23. [DOI] [PubMed] [PMC]
- 39. Wall DM, Thompson M, Mileshkin L, Loveland B, Coverdale J, Xing PX, et al. Tracking Dendritic Cells In Vivo in Patients with Multiple Myeloma. Blood. 2005;106:3469. [DOI]
- 40. Zhuang Q, Cai H, Cao Q, Li Z, Liu S, Ming Y. Tolerogenic Dendritic Cells: The Pearl of Immunotherapy in Organ Transplantation. Front Immunol. 2020;11:552988. [DOI] [PubMed] [PMC]
- 41. Cechim G, Chies JAB. *In vitro* generation of human monocyte-derived dendritic cells methodological aspects in a comprehensive review. An Acad Bras Cienc. 2019;91. [DOI]
- 42. Armitage LH, Khosravi-Maharlooei M, Meacham A, Butfiloski EJ, Viola R, Egli D, et al. High-Yield Monocyte, Macrophage, and Dendritic Cell Differentiation from Induced Pluripotent Stem Cells. bioRxiv 2021.04.29.441947 [Preprint]. 2021 [cited 2024 April 8]. Available from: https://www.bior xiv.org/content/10.1101/2021.04.29.441947v1
- 43. Cunningham S, Hackstein H. Rapid generation of monocyte-derived antigen-presenting cells with dendritic cell-like properties. Transfusion. 2021;61:1845–55. [DOI] [PubMed]
- 44. Gardner A, de Mingo Pulido Á, Ruffell B. Dendritic Cells and Their Role in Immunotherapy. Front Immunol. 2020;11:924. [DOI] [PubMed] [PMC]
- 45. Cheever MA, Higano CS. PROVENGE (Sipuleucel-T) in prostate cancer: the first FDA-approved therapeutic cancer vaccine. Clin Cancer Res. 2011;17:3520–6. [DOI] [PubMed]
- 46. Döring M, Blees H, Koller N, Tischer-Zimmermann S, Müsken M, Henrich F, et al. Modulation of TAPdependent antigen compartmentalization during human monocyte-to-DC differentiation. Blood Adv. 2019;3:839–50. [DOI] [PubMed] [PMC]
- 47. Fisher JD, Zhang W, Balmert SC, Aral AM, Acharya AP, Kulahci Y, et al. In situ recruitment of regulatory T cells promotes donor-specific tolerance in vascularized composite allotransplantation. Sci Adv. 2020;6:eaax8429. [DOI] [PubMed] [PMC]
- 48. Iwahara N, Hotta K, Iwami D, Tanabe T, Tanaka Y, Ito YM, et al. Analysis of T-cell alloantigen response *via* a direct pathway in kidney transplant recipients with donor-specific antibodies. Front Immunol. 2023;14:1164794. [DOI] [PubMed] [PMC]
- 49. Siu JHY, Surendrakumar V, Richards JA, Pettigrew GJ. T cell Allorecognition Pathways in Solid Organ Transplantation. Front Immunol. 2018;9:2548. [DOI] [PubMed] [PMC]
- Smyth LA, Lechler RI, Lombardi G. Continuous Acquisition of MHC:Peptide Complexes by Recipient Cells Contributes to the Generation of Anti-Graft CD8<sup>+</sup> T Cell Immunity. Am J Transplant. 2017;17: 60–8. [DOI] [PubMed] [PMC]

- 51. Hu C, Wu Z, Li L. Mesenchymal stromal cells promote liver regeneration through regulation of immune cells. Int J Biol Sci. 2020;16:893–903. [DOI] [PubMed] [PMC]
- 52. Pilat N, Lefsihane K, Brouard S, Kotsch K, Falk C, Steiner R, et al. T- and B-cell therapy in solid organ transplantation: current evidence and future expectations. Transpl Int. 2021;34:1594–606. [DOI] [PubMed]
- 53. Louis K, Macedo C, Bailly E, Lau L, Ramaswami B, Marrari M, et al. Coordinated Circulating T Follicular Helper and Activated B Cell Responses Underlie the Onset of Antibody-Mediated Rejection in Kidney Transplantation. J Am Soc Nephrol. 2020;31:2457–74. [DOI] [PubMed] [PMC]
- 54. Kern B, Martin F, Mengwasser J, Reutzel-Selke A, Ashraf MI, Fuehrer K, et al. 427.1: Donor Dendritic Cell Depletion Allows Novel Insights Into the Alloimmune Response Following Murine Composite Tissue and Skin Transplantation. Transplantation. 2022;106:S493–4. [DOI]
- 55. Morelli AE, Thomson AW. Tolerogenic dendritic cells and the quest for transplant tolerance. Nat Rev Immunol. 2007;7:610–21. [DOI] [PubMed]
- 56. Tanriver Y, Ratnasothy K, Bucy RP, Lombardi G, Lechler R. Targeting MHC class I monomers to dendritic cells inhibits the indirect pathway of allorecognition and the production of IgG alloantibodies leading to long-term allograft survival. J Immunol. 2010;184:1757–64. [DOI] [PubMed]
- 57. Morelli AE, Larregina AT. Concise Review: Mechanisms Behind Apoptotic Cell-Based Therapies Against Transplant Rejection and Graft versus Host Disease. Stem Cells. 2016;34:1142–50. [DOI] [PubMed] [PMC]
- 58. Moreau A, Kervella D, Bouchet-Delbos L, Braudeau C, Saïagh S, Guérif P, et al. A Phase I/IIa study of autologous tolerogenic dendritic cells immunotherapy in kidney transplant recipients. Kidney Int. 2023;103:627–37. [DOI] [PubMed]
- 59. Torres-Aguilar H, Sánchez-Torres C, Jara LJ, Blank M, Shoenfeld Y. IL-10/TGF-beta-treated dendritic cells, pulsed with insulin, specifically reduce the response to insulin of CD4+ effector/memory T cells from type 1 diabetic individuals. J Clin Immunol. 2010;30:659–68. [DOI] [PubMed]
- 60. Stax AM, Gelderman KA, Schlagwein N, Essers MC, Kamerling SWA, Woltman AM, et al. Induction of donor-specific T-cell hyporesponsiveness using dexamethasone-treated dendritic cells in two fully mismatched rat kidney transplantation models. Transplantation. 2008;86:1275–82. [DOI] [PubMed]
- 61. Flamann C, Peter K, Kreutz M, Bruns H. Regulation of the Immune Balance During Allogeneic Hematopoietic Stem Cell Transplantation by Vitamin D. Front Immunol. 2019;10:2586. [DOI] [PubMed] [PMC]
- 62. Comi M, Amodio G, Passeri L, Fortunato M, Sio FRSd, Andolfi G, et al. Generation of Powerful Human Tolerogenic Dendritic Cells by Lentiviral-Mediated IL-10 Gene Transfer. Front Immunol. 2020;11: 1260. [DOI] [PubMed] [PMC]
- Zhang M, Zheng Y, Sun Y, Li S, Chen L, Jin X, et al. Knockdown of NEAT1 induces tolerogenic phenotype in dendritic cells by inhibiting activation of NLRP3 inflammasome. Theranostics. 2019;9: 3425–42. [DOI] [PubMed] [PMC]
- 64. Wu J, Zhang H, Zheng Y, Jin X, Liu M, Li S, et al. The Long Noncoding RNA MALAT1 Induces Tolerogenic Dendritic Cells and Regulatory T Cells *via* miR155/Dendritic Cell-Specific Intercellular Adhesion Molecule-3 Grabbing Nonintegrin/IL10 Axis. Front Immunol. 2018;9:1847. [DOI] [PubMed] [PMC]
- 65. Zahorchak AF, DeRiggi ML, Muzzio JL, Sutherland V, Humar A, Lakkis FG, et al. Manufacturing and validation of Good Manufacturing Practice-compliant regulatory dendritic cells for infusion into organ transplant recipients. Cytotherapy. 2023;25:432–41. [DOI] [PubMed]
- 66. Hopewell EL, Cox C. Manufacturing Dendritic Cells for Immunotherapy: Monocyte Enrichment. Mol Ther Methods Clin Dev. 2020;16:155–60. [DOI] [PubMed] [PMC]
- 67. Li JSY, Robertson H, Trinh K, Raghubar AM, Nguyen Q, Matigian N, et al. Tolerogenic dendritic cells protect against acute kidney injury. Kidney Int. 2023;104:492–507. [DOI] [PubMed]

- Ezzelarab MB, Lu L, Shufesky WF, Morelli AE, Thomson AW. Donor-Derived Regulatory Dendritic Cell Infusion Maintains Donor-Reactive CD4<sup>+</sup>CTLA4<sup>hi</sup> T Cells in Non-Human Primate Renal Allograft Recipients Treated with CD28 Co-Stimulation Blockade. Front Immunol. 2018;9:250. [DOI] [PubMed] [PMC]
- 69. Costa D, Picascia A, Grimaldi V, Amarelli C, Petraio A, Levi A, et al. Role of HLA matching and donor specific antibody development in long-term survival, acute rejection and cardiac allograft vasculopathy. Transpl Immunol. 2024;83:102011. [DOI] [PubMed]
- 70. Zhang Q, Chen Y, Fairchild RL, Heeger PS, Valujskikh A. Lymphoid sequestration of alloreactive memory CD4 T cells promotes cardiac allograft survival. J Immunol. 2006;176:770–7. [DOI] [PubMed]
- Suchin EJ, Langmuir PB, Palmer E, Sayegh MH, Wells AD, Turka LA. Quantifying the frequency of alloreactive T cells in vivo: new answers to an old question. J Immunol. 2001;166:973–81. [DOI] [PubMed]
- 72. Nonoyama S, Hotta K, Iwahara N, Tanabe T, Hirose T, Harada S, et al. Use of Mixed Lymphocyte Reaction Assay to Evaluate Immune Tolerance before Kidney Transplantation with an Immunosuppression-Free Protocol following Hematopoietic Stem Cell Transplantation from the Same Donor. Nephron. 2023;147:621–6. [DOI] [PubMed]
- Fu J, Khosravi-Maharlooei M, Sykes M. High Throughput Human T Cell Receptor Sequencing: A New Window Into Repertoire Establishment and Alloreactivity. Front Immunol. 2021;12:777756. [DOI]
  [PubMed] [PMC]
- 74. DeWolf S, Shen Y, Sykes M. A New Window into the Human Alloresponse. Transplantation. 2016;100: 1639–49. [DOI] [PubMed] [PMC]
- 75. Liu E, Tu W, Law HK, Lau YL. Decreased yield, phenotypic expression and function of immature monocyte-derived dendritic cells in cord blood. Br J Haematol. 2001;113:240–6. [DOI] [PubMed]
- 76. Eljaafari A, Farre A, Duperrier K, Even J, Vie H, Michallet M, et al. Generation of helper and cytotoxic CD4+T cell clones specific for the minor histocompatibility antigen H-Y, after in vitro priming of human T cells by HLA-identical monocyte-derived dendritic cells. Transplantation. 2001;71: 1449–55. [DOI] [PubMed]
- 77. Dillon SR, Lewis KE, Swanson R, Evans LS, Kornacker MG, Levin SD, et al. A Dual ICOS/CD28 Antagonist ICOSL Variant Ig Domain (vIgD) Potently Suppresses Human Mixed Lymphocyte Reactions and Human/NSG Mouse Xenograft Graft vs. Host Disease (GvHD). Biol Blood Marrow Transplant. 2018;24:S187–8. [DOI]
- 78. Eljaafari-Corbin A, Mazet S, Betuel H, Rigal D, Gebuhrer L. Detection of minor histocompatibility antigen disparities by monocyte-derived dendritic cells in HLA-identical siblings. Immunol Lett. 1997;56:462–3. [DOI]
- 79. Kwon M, Choi YJ, Sa M, Park S, Shin E. Two-Round Mixed Lymphocyte Reaction for Evaluation of the Functional Activities of Anti-PD-1 and Immunomodulators. Immune Netw. 2018;18:e45. [DOI] [PubMed] [PMC]
- 80. Sharp TA, Gress RE, Sachs DH, Rosenberg SA. Assessment of mixed lymphocyte reactivity in human bone marrow cell cultures. Transplantation. 1985;40:551–6. [DOI] [PubMed]
- 81. Shimizu I, Fudaba Y, Shimizu A, Yang Y, Sykes M. Comparison of human T cell repertoire generated in xenogeneic porcine and human thymus grafts. Transplantation. 2008;86:601–10. [DOI] [PubMed] [PMC]
- Li L, Chen X, Zhang Y, Li Q, Qi C, Fei X, et al. Toll-like receptor 2 deficiency promotes the generation of alloreactive Th17 cells after cardiac transplantation in mice. Cell Immunol. 2019;338:9–20. [DOI] [PubMed]
- 83. Lindemann M. Ex vivo assessment of cellular immune function applications in patient care and clinical studies. Tissue Antigens. 2014;84:439–49. [DOI] [PubMed]

- 84. Sato T, Deiwick A, Raddatz G, Koyama K, Schlitt HJ. Interactions of allogeneic human mononuclear cells in the two-way mixed leucocyte culture (MLC): influence of cell numbers, subpopulations and cyclosporin. Clin Exp Immunol. 1999;115:301–8. [DOI] [PubMed] [PMC]
- 85. Pissas G, Eleftheriadis T. Assessment of Humoral Alloimmunity in Mixed Lymphocyte Reaction. Bio Protoc. 2019;9:e3139. [DOI] [PubMed] [PMC]
- 86. Little CJ, Kim SC, Fechner JH, Post J, Coonen J, Chlebeck P, et al. Early allogeneic immune modulation after establishment of donor hematopoietic cell-induced mixed chimerism in a nonhuman primate kidney transplant model. Front Immunol. 2024;15:1343616. [DOI] [PubMed] [PMC]
- 87. Schwarzbich MA, Gutknecht M, Guettler L, Kanz L, Salih HR, Gruenebach F, et al. Immunosuppressive Drugs Upregulate the Immunoinhibitory Receptor Osteoactivin in Monocyte-Derived Dendritic Cells. Blood. 2011;118:3232. [DOI]
- Haidinger M, Poglitsch M, Geyeregger R, Kasturi S, Zeyda M, Zlabinger GJ, et al. A versatile role of mammalian target of rapamycin in human dendritic cell function and differentiation. J Immunol. 2010;185:3919–31. [DOI] [PubMed]
- 89. Bedrosian I, Roros JG, Xu S, Nguyen HQ, Engels F, Faries MB, et al. Granulocyte-macrophage colonystimulating factor, interleukin-2, and interleukin-12 synergize with calcium ionophore to enhance dendritic cell function. J Immunother. 2000;23:311–20. [DOI] [PubMed]
- 90. Säemann MD, Kelemen P, Böhmig GA, Hörl WH, Zlabinger GJ. Hyporesponsiveness in alloreactive Tcells by NF-kappaB inhibitor-treated dendritic cells: resistance to calcineurin inhibition. Am J Transplant. 2004;4:1448–58. [DOI] [PubMed]
- 91. Min J, Liu L, Wang J, Shang CZ, Wan YL, Yu Q, et al. Morphology and immunological function of immature dendritic cells induced by sodium butyrate in human monocytes. J Clin Rehabil Tissue Eng Res. 2007;11:7485–8.
- 92. Hoeck JV, Vanhove C, Smedt SCD, Raemdonck K. Non-invasive cell-tracking methods for adoptive T cell therapies. Drug Discov Today. 2022;27:793–807. [DOI] [PubMed]
- 93. Protheroe RE, Steward CG, Mazza G, Wraith DC. Human CD4<sup>+</sup>CD25<sup>+</sup>CD127<sup>-</sup> T Cells Show Potent Dose-Dependent Inhibition of Allogeneic DC-Driven MLRs. Blood. 2006;108:5172. [DOI]
- 94. Hlavackova E, Pilatova K, Cerna D, Selingerova I, Mudry P, Mazanek P, et al. Dendritic Cell-Based Immunotherapy in Advanced Sarcoma and Neuroblastoma Pediatric Patients: Anti-cancer Treatment Preceding Monocyte Harvest Impairs the Immunostimulatory and Antigen-Presenting Behavior of DCs and Manufacturing Process Outcome. Front Oncol. 2019;9:1034. [DOI] [PubMed] [PMC]
- 95. Bellone G, Carbone A, Smirne C, Scirelli T, Buffolino A, Novarino A, et al. Cooperative induction of a tolerogenic dendritic cell phenotype by cytokines secreted by pancreatic carcinoma cells. J Immunol. 2006;177:3448–60. [DOI] [PubMed]
- 96. Hlavackova E, Pilatova K, Fedorova L, Mudry P, Mazanek P, Sterba J, et al. Pharmacotherapy preceding monocyte harvest interferes with the potency of monocyte-derived dendritic cell anticancer vaccine. Ann Oncol. 2019;30:xi13. [DOI]
- 97. Almeida CVD, Zamame JA, Romagnoli GG, Rodrigues CP, Magalhães MB, Amedei A, et al. Treatment of colon cancer cells with 5-fluorouracil can improve the effectiveness of RNA-transfected antitumor dendritic cell vaccine. Oncol Rep. 2017;38:561–8. [DOI] [PubMed]
- 98. Neves AR, Ensina LFC, Anselmo LB, Leite KRM, Buzaid AC, Câmara-Lopes LH, et al. Dendritic cells derived from metastatic cancer patients vaccinated with allogeneic dendritic cell-autologous tumor cell hybrids express more CD86 and induce higher levels of interferon-gamma in mixed lymphocyte reactions. Cancer Immunol Immunother. 2005;54:61–6. [DOI] [PubMed] [PMC]
- 99. Radice E, Bellone G, Miranda V. Enhancement of the Immunostimulatory Functions of Ex Vivo-Generated Dendritic Cells from Early-Stage Colon Cancer Patients by Consecutive Exposure to Low Doses of Sequential-Kinetic-Activated IL-4 and IL-12. A Preliminary Study. Transl Oncol. 2015;8: 327–38. [D0I] [PubMed] [PMC]

- 100. Trepiakas R, Pedersen AE, Met O, Hansen MH, Berntsen A, Svane IM. Comparison of alpha-Type-1 polarizing and standard dendritic cell cytokine cocktail for maturation of therapeutic monocyte-derived dendritic cell preparations from cancer patients. Vaccine. 2008;26:2824–32. [DOI] [PubMed]
- 101. Fang KK, Lee JB, Zhang L. Adoptive Cell Therapy for T-Cell Malignancies. Cancers (Basel). 2022;15:
  94. [DOI] [PubMed] [PMC]
- 102. Mo F, Srinivasan M, Heslop HE, Brenner MK, Mamonkin M. Rejection-Resistant Off-the-Shelf T Cells for Adoptive Cell Therapy. Biol Blood Marrow Transplant. 2019;25:S168. [DOI]
- 103. Ganguly D, Chamilos G, Lande R, Gregorio J, Meller S, Facchinetti V, et al. Self-RNA-antimicrobial peptide complexes activate human dendritic cells through TLR7 and TLR8. J Exp Med. 2009;206: 1983–94. [DOI] [PubMed] [PMC]
- 104. Torres-Aguilar H, Blank M, Kivity S, Misgav M, Luboshitz J, Pierangeli SS, et al. Tolerogenic dendritic cells inhibit antiphospholipid syndrome derived effector/memory CD4<sup>+</sup> T cell response to β2GPI. Ann Rheum Dis. 2012;71:120–8. [DOI] [PubMed]
- 105. Ju X, Silveira PA, Hsu W, Elgundi Z, Alingcastre R, Verma ND, et al. The Analysis of CD83 Expression on Human Immune Cells Identifies a Unique CD83+-Activated T Cell Population. J Immunol. 2016; 197:4613–25. [DOI] [PubMed]
- 106. Vaknin-Dembinsky A, Balashov K, Weiner HL. IL-23 is increased in dendritic cells in multiple sclerosis and down-regulation of IL-23 by antisense oligos increases dendritic cell IL-10 production. J Immunol. 2006;176:7768–74. [DOI] [PubMed]
- 107. Tai Y, Wang Q, Korner H, Zhang L, Wei W. Molecular Mechanisms of T Cells Activation by Dendritic Cells in Autoimmune Diseases. Front Pharmacol. 2018;9:642. [DOI] [PubMed] [PMC]
- 108. Andreakos E, Smith C, Monaco C, Brennan FM, Foxwell BM, Feldmann M. IκB kinase 2 but not NF-κBinducing kinase is essential for effective DC antigen presentation in the allogeneic mixed lymphocyte reaction. Blood. 2003;101:983–91. [DOI] [PubMed]
- 109. Volchenkov R, Brun JG, Jonsson R, Appel S. In vitro suppression of immune responses using monocyte-derived tolerogenic dendritic cells from patients with primary Sjögren's syndrome. Arthritis Res Ther. 2013;15:R114. [DOI] [PubMed] [PMC]
- 110. Ríos-Ríos WdJ, Sosa-Luis SA, Torres-Aguilar H. Current advances in using tolerogenic dendritic cells as a therapeutic alternative in the treatment of type 1 diabetes. World J Diabetes. 2021;12:603–15. [DOI] [PubMed] [PMC]
- 111. Qu C, Brinck-Jensen N, Zang M, Chen K. Monocyte-derived dendritic cells: targets as potent antigenpresenting cells for the design of vaccines against infectious diseases. Int J Infect Dis. 2014;19:1–5.
   [DOI] [PubMed]
- 112. Liaqat F, Kangethe RT, Pichler R, Liu B, Huber J, Wijewardana V, et al. Determination of Vaccine Immunogenicity Using Bovine Monocyte-Derived Dendritic Cells. J Vis Exp. ;2023195:e64874. [DOI] [PubMed]
- 113. Ersland K, Wüthrich M, Klein BS. Dynamic interplay among monocyte-derived, dermal, and resident lymph node dendritic cells during the generation of vaccine immunity to fungi. Cell Host Microbe. 2010;7:474–87. [DOI] [PubMed] [PMC]
- 114. Jonny J, Putranto TA, Yana ML, Sitepu EC, Irfon R, Ramadhani BP, et al. Safety and efficacy of dendritic cell vaccine for COVID-19 prevention after 1-Year follow-up: phase I and II clinical trial final result. Front Immunol. 2023;14:1122389. [DOI] [PubMed] [PMC]
- 115. Patente TA, Pinho MP, Oliveira AA, Evangelista GCM, Bergami-Santos PC, Barbuto JAM. Human Dendritic Cells: Their Heterogeneity and Clinical Application Potential in Cancer Immunotherapy. Front Immunol. 2019;9:3176. [DOI] [PubMed] [PMC]
- 116. Milhau N, Almouazen E, Bouteille S, Hellel-Bourtal I, Azzouz-Maache S, Benavides U, et al. In vitro evaluations on canine monocyte-derived dendritic cells of a nanoparticles delivery system for vaccine antigen against Echinococcus granulosus. PLoS One. 2020;15:e0229121. [DOI] [PubMed] [PMC]

- Hu Y, Su X, Liu L, Xiang Y, Yu Q, Yi S, et al. [Isolation, amplification and identification of human natural CD4<sup>+</sup>CD25<sup>+</sup> regulatory T cells in vitro]. Zhong Nan Da Xue Xue Bao Yi Xue Ban. 2014;39: 818–24. [DOI] [PubMed]
- 118. Xie L, Liu G, Liu Y, Yu Y. *In vitro* and *In vivo* CD8<sup>+</sup> T Cell Suppression Assays. Bio Protoc. 2021;11: e4020. [DOI] [PubMed] [PMC]
- 119. Zahorchak AF, Kean LS, Tokita D, Turnquist HR, Abe M, Finke J, et al. Infusion of stably immature monocyte-derived dendritic cells plus CTLA4Ig modulates alloimmune reactivity in rhesus macaques. Transplantation. 2007;84:196–206. [DOI] [PubMed]