



# Convenient estimation of oxytetracycline and polymyxin B by a novel high-performance liquid chromatography method: development and validation

Tanu Chaudhary<sup>1</sup>, Balak Das Kurmi<sup>2</sup>, Dilpreet Singh<sup>2\*</sup>

<sup>1</sup>Department of Quality Assurance, ISF College of Pharmacy, Moga, Punjab 142001, India

<sup>2</sup>Department of Pharmaceutics, ISF College of Pharmacy, Moga, Punjab 142001, India

**\*Correspondence:** Dilpreet Singh, Department of Pharmaceutics, ISF College of Pharmacy, Moga, Punjab 142001, India. [dilpreet.daman@gmail.com](mailto:dilpreet.daman@gmail.com)

**Academic Editor:** Fernando Albericio, University of KwaZulu-Natal, South Africa; University of Barcelona, Spain

**Received:** September 19, 2022 **Accepted:** November 8, 2022 **Published:** February 24, 2023

**Cite this article:** Chaudhary T, Kurmi BD, Singh D. Convenient estimation of oxytetracycline and polymyxin B by a novel high-performance liquid chromatography method: development and validation. *Explor Drug Sci.* 2023;1:6–17. <https://doi.org/10.37349/eds.2023.00002>

## Abstract


**Aim:** The aim of this research work was to develop a validated reversed-phase (RP)-high-performance liquid chromatography (HPLC) method for simultaneous estimation of oxytetracycline (OXY) and polymyxin B (PMB) in fixed-dose combination.

**Methods:** The HPLC assay method was validated on X-Bridge C<sub>18</sub> [250 mm × 4.6 mm intradermal (i.d.), 5 μm], mobile phase consisting of acetonitrile co-elution network (ACN):water containing 0.5% (v/v) orthophosphoric acid (pH 3.5) in the ratio of 80:20 respectively. The flow rate was set at 0.9 mL/min and the column was maintained at room temperature. The RP-HPLC method was validated in terms of the calibration curve (CC), linearity and range, limit of detection (LOD), and limit of quantitation (LOQ), precision, robustness, and accuracy.

**Results:** The method was found to be linear with a concentration range of 5–25 μg/mL. Precision results showed the developed method was found to be precise with a relative standard deviation [RSD (%)] value < 2. Accuracy showed acceptable recovery of prepared concentrations as per International Conference on Harmonization (ICH) guidelines. Moreover, the developed method was found to be robust and rugged, as per specified ranges. The assay of these two drugs in marketed formulation, i.e., Terramycin® Ointment showed satisfactory recovery, as per ICH guidelines. The results proved that the method can be used for the routine-based estimation of OXY and PMB.

**Conclusions:** Linear CC were obtained with a correlation coefficient ( $R^2 > 0.99$ ) with acceptable results of accuracy and precision.





**Graphical abstract.** Ultraviolet visible and HPLC method development

## Keywords

Oxytetracyclin, polymixin B, validation high-performance liquid chromatography analysis, analysis

## Introduction

The success of a pharmaceutical product relies on its qualitative analysis with a thorough validation process [1]. With the advancements in medicinal chemistry, pharmacology, and thereby the launch of new formulations, the demand for suitable testing methods has increased a lot in recent times. The acceptability of the method of analysis lies in the convenience of application and less costly approach [2]. When no conclusive approaches are available, new methodologies for evaluating novel pharmaceutical combinations are needed. Novel procedures have been developed to analyze pharmaceutical formulations with greater precision and strength. Preliminary runs are employed to optimize and validate these approaches [3]. The analytical method is an important aspect of regulatory filing in Europe and the United States, where enormous technological breakthroughs in analytical instrumentation are achieved with an emphasis on reducing analysis time without sacrificing precision and accuracy [4]. All of the new drug molecules take a too long time to have their inclusion in standard pharmacopeia books. Standards and analytical procedures for these substances and their different combinations may not be available in pharmacopeias under these circumstances [5]. As a result, new analytical procedures for these medications must be developed to assure the identification, purity, potency, and performance of drug items [6].

The validation data to support analytical procedures is collected in a planned and systematic manner by the applicant as part of the method validation process for analytical procedures. Validation is required for all analytical procedures intended for use in the analysis of clinical samples. Analytical methods are validated in accordance with International Conference on Harmonization (ICH) guidelines [7].

Oxytetracycline (OXY) has been engaged as an antimicrobial agent in the management of mastitis in dairy cows. As a consequence, (1) OXY is found in fairly large amounts in pasteurized milk along with other antibiotic residues. Similarly, honey is observed with considerable contents of OXY. OXY is used in the treatment of European foulbrood disease and American foulbrood disease due to *Paenibacillus larvae* and

*Streptococcus pluton* bacteria, respectively. (2) The honey containing these antibiotic residues is not meant for sale in advanced countries in particular [8]. (3) Further, OXY has got a wide application in dealing with fish bacterial diseases induced by *Aeromonas salmonicida*, *Aeromonas hydrophila*, and *Aeromonas sobia*. (4) Polymixin B (PMB) is also used in the treatment of various bacterial diseases associated with aquatic life. The use of antibiotics becomes unavoidable in certain circumstances, but the setting up of limits on the food contents may allow these items to be consumed safely by humans. Thus, there is an immense need for accurate and precise estimation of these antibiotics in different food samples. The capacity of an analytical process to produce test results that are directly proportional to the concentration (quantity) of analyte in the sample (within a specific range) is known as linearity. The robustness of an analytical procedure is a measure of its ability to remain unaffected by modest but deliberate changes in method parameters, and it indicates its dependability in everyday use [9]. Combination products, commonly known as fixed-dose drug combinations (FDCs), are single dose pharmaceutical formulations that combine two or more active medications. A new investigational drug in conjunction with an existing drug or two or more novel entities could also be included in the FDC formulation. FDC drugs are frequently utilized to treat single or numerous diseases. Patients are benefited from FDC products because they simplify dosing and lower prescription costs [9].

OXY is a wide-spectrum antibiotic that belongs to the tetracycline family [10, 11]. However, as these medications leave residues in milk or meat that might be directly poisonous or induce allergic reactions in some hypersensitive individuals, their use has become a severe problem in terms of infectious illnesses. Even more importantly, low-dose antibiotics in foods taken over long periods can cause problems with the proliferation of drug-resistant bacteria [12]. Different techniques are commonly used as screening methods to identify tetracycline residues in food [13], but these approaches lack specificity, give only semi-quantitative assessments of residues identified and can produce false positives. High-performance liquid chromatography (HPLC) with diode array detector (DAD) [14] and fluorimetric [15] detection has been employed for quantitative tetracycline determinations, but some chromatographic techniques are suited for multi-residue analysis [16, 17]. OXY and PMB are available as a fixed-dose topical gel formulation (0.5% strength) for cutaneous infections. The simultaneous presence of OXY and PMB in different food items and FDCs comprising these antibiotics insist on the requirement of a novel method of analysis for the efficient detection of these two drugs.

OXY and PMB have not been analyzed simultaneously till date. To reduce the analytical error and offer readers a good approach for determining the two antibiotics in a single injection, it hypothesized that these two drugs were estimated simultaneously. In this work, the authors have developed a simple, rapid, and precise method for the quantitative and qualitative measurement of both OXY and PMB in different samples.

## Materials and methods

### Chemicals and reagents

OXY (3288-A) and PMB (2332-B) pure Active Pharmaceutical Ingredients were obtained as gift samples from H.L. Pharma Private Limited (India). Acetonitrile (HPLC grade) was sourced from Loba Chemie, Mumbai, while deionized water was sourced from ISF Analytical Laboratory and also met the requirements for HPLC grade. In addition, 98% (w/v) ortho-phosphoric acid (OPA, 1772-A-212) was acquired from Loba Chemie in Mumbai. Methanol (282883211) was obtained (HPLC grade) was purchased from Rankem Chemicals Private Limited, India.

### Determination of isosbestic point by ultraviolet visible spectrophotometer

The OXY and PMB ultraviolet visible (UV) spectrophotometer (UV1900i, Shimadzu, India) were determined in Acetonitrile (ACN) by preparing a standard solution comprising 10 µg/mL of each medication. Before taking the spectrophotometer, the baseline was adjusted using the blank. Spectrophotometer was obtained three times (triplicate), and the maximum wavelength for OXY was observed to be 316 nm, while PMB's maximum wavelength was observed at 248 nm. In order to determine the isosbestic point, standard solutions

of both drugs were made with a concentration of 10 µg/mL and absorbance was noted at cross wavelengths. The isosbestic point was determined by a simultaneous equation method.

## Reversed-phase-HPLC method development

### Instrumentation and chromatographic conditions

A 515 HPLC pump and a 2489 UV/Visible spectrophotometer detector (Waters) constituted the HPLC system was utilized in the present study. C<sub>18</sub> analytical column (250 mm × 4.6 mm; particle size 5 µm; Waters) with a temperature of 42°C was used to separate the drugs. Acetonitrile and HPLC grade water in different ratios were filtered and degassed, and the pH was corrected to 3.5 using OPA and sodium hydroxide (A876721). An isocratic gradient was selected for the final analysis. For both drugs, the effluent was measured at an isosbestic point of 303 nm. Several trials were carried out to set the method varying the mobile phase composition, flow rate, and pH during method development. Empower software (Version 2.1221) was used to collect the data. The optimized chromatographic conditions for the development of the OXY and PMB HPLC method has been given in Table 1.

**Table 1.** Optimization of chromatographic conditions for reversed-phase (RP)-HPLC method

Chromatographic parameters	Variations	Peak shape
Mobile phase composition (v/v)		
ACN:Methanol:Water	50:30:20	No ideal peak
ACN:Water	50:50	Broad peak
	60:40	Not symmetric
	70:30	Sharp resolved peak, tailing
	80:20	Well resolved peak
pH variation		
ACN:Orthophosphoric acid (pH 4.0);	4.0	Broad peaks
optimized ratio (80:20)	3.0	Resolved sharp peak, tailing
	3.5	Sharp, well separated peaks

### Preparation of calibration curve

The calibration curves (CCs) of standard solutions were prepared by diluting a suitable quantity of OXY stock solution with ACN as diluent. The CC standard solutions were prepared in concentrations of 5–25 µg/mL for both PMB and OXY, respectively. Baseline correction with blank was done before taking the spectrophotometer. The absorbance was measured at  $\lambda_{\max}$  of OXY in the same solvent mixture against a blank solution using a UV spectrophotometer. The mean absorbance was measured and a graph was plotted between concentration (x-axis) and absorbance (y-axis). The linear equation and correlation coefficient were calculated for the CC.

### Method validation study

The system suitability parameter gives assurance that the analytical system is working properly and may give accurate and precise results [16]. Following optimization and development of the chromatographic method was done as per ICH Q2 (R1), validation of the analytical method has been done as per ICH Q2 (B). Different analytical tests like specificity, linearity, limit of detection (LOD), limit of quantitation (LOQ), accuracy, precision, robustness, and ruggedness as per standard guidelines [8].

### LOD and LOQ

In order to produce samples for the LOD and LOQ, the OXY and PMB stock solutions were serially diluted (as prepared). The limits of blank, quantification, and detection were determined by ACN utilizing the specified formulas and calculating them based on the standard deviation (SD) and slope of the calibration plot with these formulas 1 and 2:

$$\text{LOD} = 3.3 \times \sigma/S \quad (1)$$

$$\text{LOQ} = 10 \times \sigma/S \quad (2)$$

Where,  $\sigma$  = the SD of the response; S = slope of the CC [8].

### Linearity

Results from test samples that are linear within a specified range are directly proportional to the analyte concentration. The CC has been developed utilizing a linear progression of multiple concentrations, ranging from 5–25  $\mu\text{g/mL}$ , for both OXY and PMB, respectively. By automatically zeroing the intercept at the slope's vertices, linear regression analysis was used to determine the linearity. The least squares technique of analysis's correlation coefficient ( $R^2$ ) must not be  $< 0.990$  in order to meet the acceptance criteria.

### Intraday precision (reproducibility) and interday precision (repeatability)

By producing replicates of three distinct test concentrations (5  $\text{g/mL}$ , 15  $\text{g/mL}$ , and 25  $\text{g/mL}$ ) at 100% level, precision experiments (intraday and interday) for spectrophotometric investigations were carried out to determine the repeatability of the proposed approach. The drug amount was determined for both intraday and interday accuracy. Results were presented as relative SD [RSD (%)], with an acceptable value of  $< 2$  being regarded as statistically significant and consistent with the suggested method's capacity to be replicated.

### Ruggedness

Ruggedness is used to quantify the reproducibility of test results when laboratory test settings are changed, for as by using various analysts, different days, different reagents, etc. Three replicates at various concentrations for two separate analysts were created and analyzed for HPLC investigations. Results were reported as RSD (%) with an acceptable value of  $< 2$ , considering the procedure to be robust. The corresponding mean absorbance (UV) and peak area (HPLC) were noted.

### Robustness

Reproducibility of test findings under various temperature settings, etc., is a key component of robustness. Experiments were conducted for HPLC research by altering the flow rate, run time, and detection. The samples were made in triplicate, and the results were given as RSD (%).

### Accuracy

Accuracy is the percentage of an analyte that was added by the test from the known amount. The approach used was to prepare concentration ranges at three distinct levels (80%, 100%, and 120%) *versus* a nominal concentration (25  $\text{g/mL}$ ) for HPLC analyses. Following injection, each developed concentration was evaluated and percent recovery was calculated. In triplicate, both solutions were made and tested to calculate RSD [10, 15].

### Assay determination in marketed tablets

The method was validated to determine active content in tablet formulation. Briefly, 0.1 g of Terramycin® Ointment (Zydus Cadila Pvt Ltd, India) was taken and diluted to prepare a linear concentration of 100  $\mu\text{g/mL}$ . The prepared concentration was then analyzed by the developed UV method and simultaneously evaluated for active drug content present in the marketed formulation.

## Results

### Determination of isosbestic point by UV spectrophotometer


During method development, a diluted concentration of OXY and PMB in ACN (10  $\mu\text{g/mL}$ ) was prepared and scanned at the beer lambert range. The isosbestic point of 303 nm was observed for both drugs, respectively, using the simultaneous equation method.

### HPLC method development

#### Optimization conditions

The data of different trial approaches was given in Table S1–5. The optimal chromatographic settings for the present work are depicted in Table 1. Different trials were done to optimize the selection and concentration

of mobile phase and optimal flow rate. In the final optimization, acetonitrile and water (pH of 3.5 adjusted with OPA) were chosen, because it was shown to offer excellent peak symmetry for both of the drugs being tested, respectively. The representative chromatograms were shown in [Figure 1](#).




**Figure 1.** HPLC chromatogram. (A) OXY; (B) PMB; (C) simultaneous estimation of both OXY and PMB

#### Linearity and range

Experimental findings observed that there was a linear association in the concentration range of 5–25 µg/mL for both OXY and PMB ([Figure S1](#)). The linear concentrations developed in the investigation were well-fitted to Beer's law. In linearity curves, the correlation coefficient and equation for OXY and PMB were,  $R^2 = 0.991$  and  $R^2 = 0.996$  and the corresponding values for the y-axis were,  $y = 3,280.7x - 9,906.6$  and  $y = 2,437.6x - 10,926$ . The overlay plots of linear chromatograms of evaluated concentrations are given in [Figure 2](#).





**Figure 2.** Linearity curve chromatograms of OXY and PMB in analytical media

### LOD and LOQ

According to the final calculations, the LOD and LOQ of the current method validation process are also reported. It was observed that the LOD for OXY and PMB was 0.599 µg/mL and 0.525 µg/mL, respectively.

### Accuracy studies

For the standard guidelines, the method's sensitivity towards the analyte was found to be 100% ± 10% within the acceptability criteria, which indicates a good degree of sensitivity. The recovery study's findings are depicted in Table 2.

**Table 2.** Recovery studies of OXY and PMB by HPLC method

Drugs	Nominal concentration (µg/mL)	Level of addition (%)	Concentration prepared (µg/mL)	Amount recovered (µg/mL)	RSD (%)	Recovery (%)
OXY	10	80	8	7.55	1.63	94.41
	10	100	10	9.44	1.27	95.6
	10	120	12	11.53	1.14	97.25
PMB	10	80	8	7.44	1.67	94.62
	10	100	10	9.47	1.34	96.2
	10	120	12	11.7	1.27	99.08

### Intraday precision (reproducibility) and interday precision (repeatability)

Three replicates of 100 percent test concentrations of OXY and PMB, ranging from 5–25 µg/mL were evaluated and results are reported in Tables 3–4 and Tables 5–6. The developed concentration shows precise results as all the evaluated concentration shows < 2 RSD (%) over a recommended validation parameter. Results concluded that the developed method is found to be precise in nature.

**Table 3.** Interday precision and method precision studies of OXY

Drug	Concentration (µg/mL)	Day 1		Day 2		Day 3	
		Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)
OXY	5	10,157	1.442	10,078	1.7	9,866	1.214
	10	20,153	1.437	20,092	1.929	19,139	1.909
	15	40,054	1.241	39,195	1.677	39,493	1.587
	20	55,194	1.577	55,360	1.509	54,994	1.507
	25	71,819	1.302	71,259	1.313	71,028	1.141
Mean	-	-	1.399	-	1.625	-	1.471
Overall RSD (%)		1.498					

-: blank cell

**Table 4.** Intraday precision and method precision studies of OXY

Drug	Concentration (µg/mL)	Morning		Afternoon		Evening	
		Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)
OXY	5	10,009	1.285	9,984	1.205	9,997	1.185
	10	19,925	1.175	19,812	0.949	20,076	0.945
	15	39,080	1.321	39,135.7	1.103	39,108	1.026
	20	54,887	0.511	55,263.3	0.683	54,933	0.574
	25	71,062	0.597	71,192.7	0.530	71,113	0.471
Mean	-	-	0.977	-	0.894	-	0.84
Overall RSD (%)							0.903

-: blank cell

**Table 5.** Interday precision and method precision studies of PMB

Drug	Concentration (µg/mL)	Day 1		Day 2		Day 3	
		Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)
PMB	5	2,892	0.856	2,921	0.977	2,864	0.918
	10	12,027	0.838	11,850	0.843	11,800	0.856
	15	24,487	0.981	24,649	0.922	24,599	0.991
	20	39,297	0.890	39,504	0.878	39,512	0.870
	25	50,486	0.711	50,507	0.863	50,328	0.970
Mean	-	-	0.855	-	0.896	-	0.921
Overall RSD (%)							0.890

-: blank cell

**Table 6.** Intraday precision and method precision studies of PMB

Drug	Concentration (µg/mL)	Morning		Afternoon		Evening	
		Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)	Peak area (n = 3)	RSD (%)
PMB	5	2,864	0.918	2,893	0.909	2,854	0.957
	10	11,800	0.856	11,877	0.875	11,879	0.88
	15	24,599	0.991	24,649	0.922	24,606	0.939
	20	39,512	0.87	39,504	0.878	39,567	0.997
	25	50,328	0.97	50,507	0.863	50,437	0.751
Mean	-	-	0.921	-	0.889	-	0.904
Overall RSD (%)							0.904

-: blank cell

## Robustness

There are a number of ways to determine robustness, including varying flow rate ( $\pm 0.2$  mL/min) and adjusting pH (3.5–4.0), adjusting wavelength ( $\pm 2$  nm), and evaluating how this affects method sensitivity and efficacy. The final RSD (%) values of robustness studies after changing the minor proportions in optimized conditions are given in [Tables 7–8](#). Results indicated that the RSD (%) values were in the acceptable ranges, suggesting that the developed method was found to be rugged [18].

**Table 7.** Summary of OXY and PMB results by robustness method

Parameters	Attributes		301 nm	303 nm	305 nm	Mean	± SD	RSD (%)
Variable conditions for OXY	Effect of change in wavelength	5 µg/mL	10,279	9,989	9,984	10,084	168.893	1.674
		15 µg/mL	38,674	39,998	39,554	39,408.67	673.858	1.709
		25 µg/mL	70,468	71,885	71,926	71,426.33	830.194	1.162
	Mean							1.515
Variable conditions for PMB	Effect of change in wavelength	5 µg/mL	2,871	2,911	2,857	2,879.67	28.023	0.973
		15 µg/mL	24,313	24,283	24,656	24,417.3	207.235	0.848
		25 µg/mL	50,889	50,225	50,985	50,699.7	413.866	0.816
	Mean							0.879



**Table 8.** Summary of OXY and PMB results by ruggedness method

Parameters	Attributes		0.7 (mL/min)	0.9 (mL/min)	1.1 (mL/min)	Mean	± SD	RSD (%)
Variable conditions for OXY	Effect of change in flow rate (mL/min)	5 µg/mL	9,969	9,899	10,224	10,030.67	171.05	1.705
		15 µg/mL	38,504	38,610	39,754	38,956	693.117	1.779
		25 µg/mL	70,268	71,991	71,984	71,414.33	992.76	1.39
	Mean							1.624
	Overall RSD (%)							1.569
Variable conditions for PMB	Effect of change in flow rate (mL/min)	5 µg/mL	2,825	2,846	28,770	2,849.33	26.159	0.918
		15 µg/mL	24,323	24,384	24,743	24,483.33	226.936	0.926
		25 µg/mL	50,965	50,195	50,376	50,512	402.612	0.797
	Mean							0.880
	Overall RSD (%)							0.879

### Assay determination

For the marketed formulation, an assay was performed to check the purity of each drug in the formulation and the percentage purity of the drug was calculated. Percentage purity for OXY and PMB was found to be 99.6% and 89%, respectively. On the consideration of the value of SD, there is no interference of excipients in the assessment of drugs [19]. The summary of validation parameters is given in Table 9. The results of the assay determination are shown in Table 10.

**Table 9.** Summary of validation parameters for HPLC method

Validation parameters	OXY	PMB
Absorption maxima ( $\lambda_{\max}$ ; nm)	303	
Linearity and range (µg/mL)	5–25	
Coefficient of determination ( $R^2$ )	0.9915	0.9961
Regression equation (y)	3,280.7x – 9,906.6	2,437.6x – 10,926
Slope (m)	3,280.7	2,437.6
Intercept (c)	9,906	10,926
LOD (µg/mL)	0.197	0.173
LOQ (µg/mL)	0.599	0.525
Precision	1.21	0.89
Robustness	1.569	0.879

**Table 10.** Assay of OXY and PMB by HPLC method

Drug name	Drug claim (mg)	Amount found (mg)	Estimation (%)
OXY	5	4.98	99.6
PMB	1	0.89	89

## Discussion

The liquid chromatographic method was designed and fine-tuned to ensure reproducibility and accuracy. Both the hit-and-trial approach and the polarity of the mobile phase were used as selection criteria in the process. In the early stages of the measurement of OXY and PMB, a number of solvents in varying proportions were utilized [17, 20]. Although several mixtures of organic solvents, either used in 100% concentration or in conjunction with water in varying proportions, were tried, however, the separation did not meet the required mark. Following system suitability parameters were considered for mobile phase selection which is written below:

- (1) A sharp peak was observed
- (2) Tailing and fronting were not observed
- (3) Appropriate elution time
- (4) Suitable number of theoretical plates and home-exercise therapy program (HETP)

As the longitudinal broadening is inversely proportional to the flow rate of the mobile phase system, maintaining the uniformity of the flow rate is highly important [21]. When the flow rate is too high or too low, it has an effect on the Gaussian peak and can lead to symmetry problems in the overall peak. Hence, a low mobile phase composition is able to reduce peak broadening and tailing and able to improve theoretical plates. During the final analysis, ACN: water (pH 3.5 adjusted with orthophosphoric acid) in the ratio of 80:20 (v/v) was the selected and optimal mobile phase, as it gives excellent separation of molecules and detection [22].

For a particular range of analyte concentration, a procedure's linearity refers to its capacity to produce test findings that are directly proportional to the sample's concentration [23]. According to the least squares method of analysis, the developed method was shown to be linear. LOD is the ability of the method to detect the lowest amount of drug and LOQ is the ability of the method to quantify the lowest amount of drug [24]. According to guidelines, LOD and LOQ are observed to be within uniform limits. Recovery experiments were conducted at 80%, 100%, and 120% of the test concentrations of both drugs to verify the method's accuracy. It was necessary to produce the samples in triplicate and to compute the percent recovery for each sample [25]. Moreover, the robustness and ruggedness define the toughness and high sensitivity of the developed method as it doesn't show any fluctuations through minuscule fluctuations [25, 26]. Analysis revealed RSD (%) values < 2% for the performed studies which define the corrected validation of the developed method as per ICH guidelines. For the assay determination, the active content of the marketed formulation showed an acceptable range of 90–110% in the ointment which successfully tested the acceptability of the developed method for routine analysis.

In conclusion, the RP-HPLC method has been used for the estimation of OXY and PMB in a bulk material. The method was finally optimized as: mobile Phase; ACN:water (80:20) with pH 3.5 using 0.5% (v/v) orthophosphoric acid, and flow rate: 0.9 mL/min. For both OXY and PMB, validation studies were done and RSD (%) was found to be < 2%. For the marketed product, an assay was performed to check the purity of OXY and PMB in the marketed product and the percentage purity of both drugs was calculated. Percentage estimation for OXY was found to be 98.64%, and for PMB it was found to be 90.84%. The developed HPLC method for determining OXY and PMB was simple, specific, accurate, precise, rapid, and economical, indicating its adequacy for routine pharmaceutical analysis of OXY and PMB.

## Abbreviations

ACN: acetonitrile

CC: calibration curve

FDCs: fixed-dose drug combinations

HPLC: high-performance liquid chromatography

ICH: International Conference on Harmonization

LOD: limit of detection

LOQ: limit of quantitation

OPA: ortho-phosphoric acid

OXY: oxytetracycline

PMB: polymyxin B

RP: reversed-phase

RSD: relative standard deviation

SD: standard deviation

UV: ultraviolet visible

## Supplementary materials

The supplementary material for this article is available at: [https://www.explorationpub.com/uploads/Article/file/10082\\_sup\\_1.pdf](https://www.explorationpub.com/uploads/Article/file/10082_sup_1.pdf).

## Declarations

### Author contributions

DS: Conceptualization, Investigation, Writing—original draft. TC: Conceptualization, Investigation, Writing—original draft. BDK: Validation, Writing—review & editing, Supervision.

### Conflicts of interest

The authors declared that they have no conflict of interest.

### Ethical approval

Not applicable.

### Consent to participate

Not applicable.

### Consent to publication

Not applicable.

### Availability of data and materials

Not applicable.

### Funding

Not applicable.

### Copyright

© The Author(s) 2023.

## References

1. Budhavant SU, Shinde KC, Sable AR, Pawar SD. A review: introduction and over view study on pharmaceuical process of validation. *World J Pharm Res.* 2020;9:640–62.
2. Thapar S, Chawla A, Gupta GK. A concise review on method development and validation parameters. *GSC Biol Pharm Sci.* 2022;19:66–77.
3. Willard HH, Merritt LL Jr, Dean JA, Settle FA Jr. Instrumental methods of analysis, 7th edition. United States: Wadsworth Publishing Company; 1988.
4. Matuszewski BK, Constanzer ML, Chavez-Eng CM. Strategies for the assessment of matrix effect in quantitative bioanalytical methods based on HPLC – MS/MS. *Anal Chem.* 2003;75:3019–30.
5. Burgess C. Valid analytical methods and procedures: a best practice approach to method selection. Royal Society of Chemistry; 2000.
6. Ermer J. Validation in pharmaceutical analysis. Part I: an integrated approach. *J Pharm Biomed Anal.* 2001;24:755–67.
7. Huber L, George SA. Diode array detection in HPLC. *Chemistry*; 1993.
8. ich.org [Internet]. Geneva: The International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use; c2022 [cited 2022 May 3]. Available from: <https://www.ich.org>
9. Gautam CS, Saha L. Fixed dose drug combinations (FDCs): rational or irrational: a view point. *Br J Clin Pharmacol.* 2008;65:795–6.
10. Fletouris DJ, Psomas JE, Botsoglou NA. Trace analysis of oxytetracycline and tetracycline in milk by high-performance liquid chromatography. *J Agric Food Chem.* 1990;38:1913–7.
11. Dihuidi K, Kucharski MJ, Roets E, Hoogmartens J, Vanderhaeghe H. Quantitative analysis of doxycycline and related substances by high-performance liquid chromatography. *J Chromatogr A.* 1985;325:413–24.

12. Levy SB. Antibiotic use for growth promotion in animals: ecologic and public health consequences. *J Food Prot.* 1987;50:616–20.
13. Kurittu J, Lönnberg S, Virta M, Karp M. Qualitative detection of tetracycline residues in milk with a luminescence-based microbial method: the effect of milk composition and assay performance in relation to an immunoassay and a microbial inhibition assay. *J Food Prot.* 2000;63:953–7.
14. Furusawa N. Rapid liquid chromatographic determination of oxytetracycline in milk. *J Chromatogr A.* 1999;839:247–51.
15. De Wasch K, Okerman L, De Brabander H, Van Hoof J, Croubels S, De Backer P. Detection of residues of tetracycline antibiotics in pork and chicken meat: correlation between results of screening and confirmatory tests. *Analyst.* 1998;123:2737–41.
16. Gamal M, Naguib IA, Abdelfatah RM. Simultaneous analysis of oxytetracycline hydrochloride, lidocaine, and bromhexine hydrochloride in the presence of many interfering excipients. *Arch Pharm (Weinheim).* 2021;354:e2100131.
17. Sokol J, Matisova E. Determination of tetracycline antibiotics in animal tissues of food-producing animals by high-performance liquid chromatography using solid-phase extraction. *J Chromatogr A.* 1994;669:75–80.
18. Temboot P, Kaewpaiboon S, Tinpun K, Nakpeng T, Khalil R, Ul-Haq Z, et al. Potential of sodium deoxycholate sulfate as a carrier for polymyxin B: physicochemical properties, bioactivity and *in vitro* safety. *J Drug Deliv Sci Technol.* 2020;58:101779.
19. Cao G, Ali FE, Chiu F, Zavascki AP, Nation RL, Li J. Development and validation of a reversed-phase high-performance liquid chromatography assay for polymyxin B in human plasma. *J Antimicrob Chemother.* 2008;62:1009–14.
20. Oka H, Ikai Y, Ito Y, Hayakawa J, Harada KI, Suzuki M, et al. Improvement of chemical analysis of antibiotics. XXIII. Identification of residual tetracyclines in bovine tissues by electrospray high-performance liquid chromatography–tandem mass spectrometry. *J Chromatogr B Biomed Sci Appl.* 1997;693:337–44.
21. Cinquina AL, Longo F, Anastasi G, Giannetti L, Cozzani R. Validation of a high-performance liquid chromatography method for the determination of oxytetracycline, tetracycline, chlortetracycline and doxycycline in bovine milk and muscle. *J Chromatogr A.* 2003;987:227–33.
22. Hede K. Antibiotic resistance: an infectious arms race. *Nature.* 2014;509:S2–3.
23. Kaale E, Chambuso M, Kitwala J. Analysis of residual oxytetracycline in fresh milk using polymer reversed-phase column. *Food Chem.* 2008;107:1289–93.
24. Falagas ME, Kasiakou SK. Toxicity of polymyxins: a systematic review of the evidence from old and recent studies. *Crit Care.* 2006;10:R27.
25. Carr C Jr, Morrison DC. Mechanism of polymyxin B-mediated lysis of lipopolysaccharide-treated erythrocytes. *Infect Immun.* 1985;49:84–9.
26. Gai Z, Samodelov SL, Kullak-Ublick GA, Visentin M. Molecular mechanisms of colistin-induced nephrotoxicity. *Molecules.* 2019;24:653.